



**Allium**  
Urological Solutions

# Allium URS Covers it all



A leading minimally invasive long term solution for malignant & benign ureteric strictures, ranging from simple stricture indications to complications management such as ureteral fistulas, ureteral bridging and anastomotic strictures.



How many stents do you see in the X-Ray?




**Self-expandable**



# The URS - Ureteral Stent

New CE approved indication, Ureteric Fistula Treatment



Placement of a fully polymeric covered Nitinol stent to minimize and close the leakage of a ureteric fistula

## Selected Publications

### Multicenter Experience

“The use of Allium ureter stents is feasible, safe and effective”

### Emerging Market-China Experience: Allium ureteral stent for ureteral stenosis

“Allium ureteral stent is an effective way... especially in patients receiving radiotherapy.”

### Product Range

Item#	Stent Model	Diameter	Length
URS-O-R-8-60	<b>Retrograde delivery system 80 cm Stent without anchor</b>	8 mm	60 mm
URS-O-R-8-80		8 mm	80 mm
URS-O-R-8-100		8 mm	100 mm
URS-O-R-8-120		8 mm	120 mm
URS-O-R-9-200		9 mm	200 mm
URS-O-R-10-60		10 mm	60 mm
URS-O-R-10-80		10 mm	80 mm
URS-O-R-10-100		10 mm	100 mm
URS-O-R-10-120		10 mm	120 mm
URS-A-8-100	<b>Antegrade delivery system 80 cm Stent with anchor</b>	8 mm	100 mm
URS-A-8-120		8 mm	120 mm
URS-A-10-100		10 mm	100 mm
URS-A-10-120		10 mm	120 mm
URS-R-8-100	<b>Retrograde delivery system 80 cm Stent with anchor</b>	8 mm	100 mm
URS-R-8-120		8 mm	120 mm
URS-R-10-100		10 mm	100 mm
URS-R-10-120		10 mm	120 mm



Scan to register to Allium Partner Zone for downloading URS publications, webinars & manuals