

# PRODUCT CATALOG



**AMSINO**  
Advancing Healthcare Worldwide®





<b>AMSure® Enteral Feeding Systems..... 1</b>	Pediatric Urine Collector.....37
PUGGLE® Enteral Feeding Pump.....2	Irrigation Trays and Kits .....38
Sentinel™ Enteral Feeding Pump.....3	Syringes.....39
ENFit® and Non-ENFit Sets.....3	Specimen Containers .....40
ENFit and Non-ENFit Syringes.....3	
ENFit and Non-ENFit Irrigation Kits and Trays.....5	
<b>AMSafe® Medication Delivery Systems ..... 6</b>	<b>AMSure Respiratory Therapy Systems.....41</b>
Gravity Administration Sets.....7	Respiratory Masks and Nasal Cannulas .....42
10 Drops/mL.....7	Nebulizer Kits.....45
15 Drops/mL.....9	Accessories.....46
60 Drops/mL.....13	Aerosol Drainage Bag.....46
AMSafe-3® Selectable Drip Chamber	Heat and Moisture Exchangers .....47
10, 15 & 60 drops/mL.....14	Tracheostomy Clean and Care Trays.....47
Secondary Piggyback Sets.....15	Sterile Water.....48
Y-type Blood Sets.....15	Suction Catheters, Trays and Kits.....49
IV Extension Sets.....16	Suction Tubing with Connectors .....51
Standard Extension Sets.....16	Suction Yankauers and Kits .....52
Microbore Extension Sets.....17	
Smallbore Extension Sets.....17	<b>Sterile Water &amp; Saline Products..... 53</b>
Flow Controller Extension Sets.....17	Sterile Water for Inhalation, USP .....54
Pressure Infusion Extension Sets.....18	Sterile Water Prefilled Humidifiers, USP .....54
Needle-Free Connectors .....18	Nebulizer Kits.....45
Bag Spike.....18	Saline Products .....55
Needle-Free Vial Adapter.....18	PreFilled Saline Flush Syringes .....55
IV Pump Administration Sets.....19	Sodium Chloride Inhalation Solution, (0.9%) USP .....56
Sure-Lok® Needle-Free IV Pump Administration Sets.....19	Device and Wound Irrigation .....56
Pre-Pierced Y-Injection IV Pump Administration Sets.....20	
<b>AMSure Urological Care ..... 21</b>	<b>AMSure Surgical &amp; Nursing Systems.....57</b>
Foley Catheter Trays & Collection.....22	Suction Tubing with Connectors .....58
Pre-Connected Foley Catheter Urine Drainage Bag Trays...22	Suction Yankauer and Kits .....59
Add-a-Foley Catheter Trays.....23	Suction Device Holster and Kits.....59
Pre-Connected Foley Catheter Urine Meter Kits .....24	Specimen Containers.....60
Add-a-Foley Urine Meter Trays.....25	Vaginal Speculums.....60
Foley Insertion Trays.....25	Disposable Umbilical Cord Clamp.....60
Urinary Drainage Bags .....26	Ear/Ulcer Syringes.....61
Urine Meter.....27	Cleansing Enema Sets.....61
ACCU-MEASURE® Pediatric Urine Meter.....27	Connectors.....62
Foley Catheters.....28	
Hydrophilic Latex Foley Catheters.....28	<b>Suction &amp; Waste Fluid Management Systems ..... 63</b>
Silicone-Coated Latex Foley Catheters.....29	RECEPTAL® Outer Shells .....64
100% Silicone Foley Catheters.....30	RECEPTAL Liners.....65
Intermittent Urinary Catheters .....31	RECEPTAL Twist Liners.....66
Urethral Insertion Trays and Kits .....31	iRECEPTAL® Surgical Waste Fluid Management System .....67
Closed System Urethral Drainage Kits.....32	EZE-VAC® Suction Canister System .....69
Self-Catheterization Kits with R-Polished™ Eye Catheters...33	EZE-VAC 2.0 Suction Canister System.....70
PVC Intermittent Urethral Catheter with R-Polished Eyes...34	ACCU-MEASURE Critical Measurement Devices.....70
PVC Intermittent Urethral Catheter .....35	RECEPTAL Solidifiers .....71
Red Rubber Latex Urethral Catheter.....35	Suction Tubing with Connectors .....72
Irrigation Sets.....36	Suction Yankauers and Kits .....73
Cystoscopy & TUR Irrigation Sets.....36	Suction Device Holster System .....74
Ancillary & Patient Care Products.....37	Accessories.....74
Urinary Leg Bags.....37	Suction Catheters and Trays.....77
Extension Tube for Leg Bags .....37	Suction Catheter Kits.....78
	Reusable Hardware.....79
	Regulator Mounts .....82
	<b>Index of REFs ..... 83</b>





## Optimizing Nutritional Delivery

AMSure Enteral Feeding Systems enhance delivery from nutrition source to patient feeding tube. The comprehensive enteral product portfolio is available with standard and ENFit® connections. The systems are designed to provide best patient outcomes and includes:

- PUGGLE® Enteral Feeding Pump .....2
- Sentinel™ Enteral Feeding Pump.....3
- ENFit & Non-ENFit Enteral Feeding Sets.....3
- ENFit & Non-ENFit Enteral Syringes .....4
- ENFit & Non-ENFit Irrigation Kits & Trays .....5

## PUGGLE® Enteral Feeding Pump

The versatile PUGGLE Enteral Feeding Pump is designed for both ambulatory and stationary applications for pediatric to adult patients. The simplified settings make feeding, whether stationary or on-the-go, a snap.

**Puggle is only available the US.**

- Large color LED screen for easy at-a-glance view of clock and settings
- Compact lightweight design fits into backpacks and transports easily
- Quiet motor with continuous peristaltic rotary pump action
- 20 hour battery life that fully charges in only 2.5 hours

**Dimensions:** 5.4 X 5.0 X 2.5 inches (13.7 x 12.7 x 6.4 cm)  
**Weight:** 1.43 lbs (649 g)  
**Operating Range:** 1 to 400 mL/h in 1 mL increments  
**Delivered Volume:** 1 – 9999 mL in 1 mL increments  
**Accuracy:** +/- 5%



REF	DESCRIPTION	QTY/UNIT
PG001	PUGGLE Enteral Feeding Pump with Pole Clamp	1/Case
PE90ST	PUGGLE Enteral Pump Spike Set, Anti-free Flow Valve, ENFit® and Transition Connector	30/Case
PE90B5T	PUGGLE Enteral Pump Set with 500 mL Bag, Anti-free Flow Valve, ENFit and Transition Connector	30/Case
PE90B12T	PUGGLE Enteral Pump Set with 1200 mL Bag, Anti-free Flow Valve, ENFit and Transition Connector	30/Case



## Sentinel™ Enteral Feeding Pump

The Sentinel enteral feeding pump is easy to use and program which minimizes training time and helps to improve programming efficacy. The AMSure® enteral feeding pump sets are recommended to be used with the Sentinel pump.

- Large display and buttons promote ease of programming
- Intuitive design simplifies tubing set-up
- Unique drop sensor minimizes nuisance alarms
- Tubing securements help support uninterrupted feeding
- Polycarbonate pump housing increases durability
- Designed for feeding or fluid applications

REF	DESCRIPTION	QTY/UNIT
S-1000-SI	Sentinel Enteral Pump	1/Each



ENF1200K

## ENFit® and Non-ENFit Enteral Feeding Sets

AMSure Enteral Feeding Sets are available with and without ENFit connectors. The feeding sets are intended for administration of enteral nutrition, fluids or medication. AMSure enteral feeding pump sets are recommended to use with the Sentinel™ enteral pump and are compatible with the Compat™ and EntraFlo™ pumps.

- ENFit connector sets are compliant with the ISO 80369-3 to help reduce risk of misconnections
- Patient label provides convenient patient documentation
- ENFit sets include pre-attached stepped transition connector
- Non-sterile, BPA-free, not manufactured with natural rubber latex

REF	ENFit CONNECTOR	FEEDING BAG SIZE			FEEDING SET TYPE		PUMP COMPATIBILITY		QTY/UNIT
		Spike Only- No Bag	500 mL	1200 mL	Gravity Set	Pump Set	Sentinel®	Compat™/ EntraFlo™	
E-1004		•				•	•	•	30/Case
E-1200				•		•	•	•	30/Case
E-1203				•	•				30/Case
ENF1004K ©	•	•				•	•	•	30/Case
ENF1200K ©	•			•		•	•	•	30/Case
ENF1203K ©	•			•	•				30/Case
ENF500K ©	•		•			•	•	•	30/Case

EntraFlo™ is a registered trademark of Medline Industries, Inc.  
 Compat™ is a registered trademark of Société des Produits Nestlé, S.A Corporation.  
 Sentinel™ is a registered trademark of Alcor Scientific Inc. ENFit® is a registered trademark of GEDSA.

## ENFit® and Non-ENFit Enteral Syringes

AMSure® Enteral Syringes are available with or without a transition adapter ENFit tip. The ENFit adapter syringe tip is compliant with the ENFit ISO 80369-3 design standard which supports the goal to improve patient safety by reducing the risk of misconnections. All syringes include a tip protector and are DEHP-free, BPA-free and not manufactured with natural rubber latex.

### Bulb Syringe


- Flexible soft bulb promotes one-handed administration

### Thumb Controlled Piston Syringe

- Designed for maximum control and one-handed administration
- Precision crafted O-ring plunger ensures a positive seal
- Re-sealable pole bag designed for safe and convenient storage

### Flat Top Piston Syringe

- Visible dosing marks designed for ease of use
- Precision crafted O-ring plunger ensures a positive seal
- Visible dosing marks designed for ease of use



REF	ENFIT TIP	TIP ADAPTER	60 ML SYRINGE TYPE			PACKAGE TYPE			BULK	STERILE	QTY/UNIT
			Bulb	Thumb Control Ring	Flat Top	Poly Pouch	Form Fill Seal	Resealable Pole Bag			
AS010			•						•		50/Case
AS011 ©			•							•	50/Case
AS115 ©		•			•	•				•	50/Case
AS116 ©		•			•			•	•		30/Case
ENS015 ©	•			•		•				•	50/Case
ENS016	•			•				•	•		30/Case
ENS115 ©	•				•	•				•	50/Case
ENS116	•				•			•	•		30/Case
ENS115NS	•				•			•	•		30/Case
ENS115T	•*				•			•	•		30/Case

\*Includes Transition Adapter



# ENFit® and Non-ENFit Enteral Feeding Irrigation Kits and Trays

AMSure® Enteral Feeding Irrigation Kits and Trays are designed for simplicity and efficiency of care.

Kits are non-sterile and include a patient label, 500 mL graduated container, either a poly pouch or re-sealable bag and a 60 mL syringe with or without an ENFit tip.

Trays are sterile and include 500 mL graduated container, 1000 mL tray, alcohol prep pad and moisture-proof under pad, a 60 mL syringe with or without an ENFit tip.

- ENFit® connector sets are compliant with ISO 80369-3 to help reduce risk of misconnections
- DEHP-free, BPA-free, not manufactured with natural rubber latex
- Includes syringe tip protector cap



## Enteral Feeding Irrigation Kits

REF	ENFIT® SYRINGE TIP	500 mL GRADUATED CONTAINER	60 mL BULB SYRINGE	60 mL THUMB CONTROL PISTON SYRINGE	POLY POUCH	RE-SEALABLE POLE BAG	PATIENT LABEL	QTY/UNIT
AS121		•	•				•	30/Case
AS126		•		•			•	30/Case
AS127		•		•		•	•	30/Case
ENK127 ©	•	•		•		•	•	30/Case

## Enteral Feeding Irrigation Trays

REF	ENFIT® SYRINGE TIP	1000 mL OUTER TRAY	500 mL GRADUATED CONTAINER	60 mL BULB SYRINGE	60 mL THUMB CONTROL PISTON SYRINGE	STERILE	ALCOHOL PAD	LARGE UNDERPAD	QTY/UNIT
AS130 ©		•	•	•		•	•	•	20/Case
AS136 ©		•	•		•	•	•	•	20/Case
ENK136 ©	•	•	•		•	•	•	•	20/Case

# AMSafe®

## MEDICATION DELIVERY SYSTEMS



## Safety First

The AMSafe IV product portfolio provides solutions to help improve patient and caregiver safety and infection prevention. Our IV product selection includes:

- Gravity Administration IV Sets.....7
  - 10 Drops/mL.....7
  - 15 Drops/mL.....9
  - 60 Drops/mL.....13
  - AMSafe-3® Selectable Drip Chamber, 10, 15 & 60 drops/mL.....14
  - Secondary Piggyback Sets .....15
  - Y-type Blood Sets .....15
- IV Extension Sets .....16
  - Standard Extension Sets.....16
  - Microbore Extension Sets.....17
  - Smallbore Extension Sets.....17
  - Flow Controller Extension Sets.....17
  - Pressure Infusion Extension Sets .....18
- Needle-Free Connectors .....18
- Bag Spike.....18
- Needle-Free Vial Adapter .....18
- IV Pump Administration Sets .....19
  - Sure-Lok® Needle-Free IV Pump Administration Sets.....19
  - Pre-Pierced Y-Injection IV Pump Administration Sets.....20



# Gravity Administration IV Sets

Amsino's comprehensive line of gravity administration IV sets are designed to improve safety and infection prevention for both caregivers and patients.

- Available in needle-free, pre-pierced and injection Y-sites in a variety of configurations to meet your clinical needs
- Available with or without pre-attached extension sets
- 10, 15 or 60 drops/mL sets available
- Includes a rotating male luer lock
- Sterile, non-pyrogenic, DEHP-free, not manufactured with natural rubber latex

## 10 Drops/mL General Primary Sets

REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
107201	10	72 in. (183 cm)	Injection Y-site	14 mL	50/Case
108301 ©	10	83 in. (211 cm)	Injection Y-site	15.6 mL	50/Case
108303	10	83 in. (211 cm)	Sure-Lok® Needle-free Y-site	16 mL	50/Case
108304	10	89 in. (226 cm)	Non-vented, Roller Clamp (2), Sure-Lok Needle-free Y-site with Pre-attached Ext Set with Pinch Clamp	17.5 mL	50/Case
109602 ©	10	96 in. (244 cm)	(2) Injection Y-sites	17.5 mL	50/Case
108306 ©	10	83 in. (211 cm)	Pre-pierced Y-site, Sure-Lok Needle-free Y-site	16 mL	50/Case
MRSA10-84	10	84 in. (213 cm)	Vented/Non-vented, Pre-pierced Y-site, MicroClave® Clear Needle-free Y-site	10.5 mL	50/Case
AC3110 ©	10	105 in. (266 cm)	(2) Sure-Lok Needle-free Y-sites	20 mL	50/Case
MRA10E-95	10	95 in. (241 cm)	Vented/Non-vented, Pre-pierced Y-site, Slide Clamp, Pre-attached Ext Set, Slide Clamp, MicroClave Clear Needle-free Y-site	18.5 mL	50/Case

## 10 Drops/mL General Primary Sets (continued)

REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
108305 ©	10	89 in. (226 cm)	Pre-pierced Y-site, Sure-Lok® Needle-free Y-site with Pre-attached Ext Set, Pinch Clamp	17.5 mL	50/Case
MRA10E-104	10	104 in. (264 cm)	Vented/Non-vented, MicroClave® Clear Needle-free Y-sites, Pre-pierced Y-site, with Pre-attached Ext Set, Slide Clamp, Pre-pierced Y-site, MicroClave Clear Needle-free Y-site	13.4 mL	50/Case
MRSA10E-107	10	107 in. (272 cm)	Vented/Non-vented, (2) MicroClave Clear Needle-free Y-sites, Pre-attached Ext Set	13.4 mL	50/Case
108308 ©	10	83 in. (211 cm)	(2) Sure-Lok Needle-free Y-sites, Pre-pierced Y-site	16.6 mL	50/Case

## 10 Drops/mL Flow Controller Sets

REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
AFS102 ©	10	100 in. (254 cm)	Vented/Non-vented, Flow Controller, Slide Clamp, Sure-Lok Needle-free Y-site	18.5 mL	50/Case
AFS208 ©	10	100 in. (254 cm)	Vented/Non-vented, Flow Controller, Slide Clamp, (2) AMSafe® Needle-free Y-sites	19.2 mL	50/Case
AFS208B	10	96 in. (244 cm)	Vented/Non-vented, Back Check Valve, Flow Controller, (2) AMSafe Needle-free Y-sites, Slide Clamp	19 mL	50/Case








# Gravity Administration IV Sets

## 10 Drops/mL



### General Primary Sets with Back Check Valve

- Provides the flexibility to infuse more than one bag of medication or solution
- Featuring a back check valve to prevent fluid exchange or backflow into the primary infusion line

REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
1010001B ©	10	100 in. (254 cm)	Vented/Non-vented, Back Check Valve, Sure-Lok® Needle-free Y-site	19.3 mL	50/Case
					
108302 ©	10	83 in. (211 cm)	Vented/Non-vented, Back Check Valve, Pre-pierced Y-site, Sure-Lok Needle-free Y-site, Rotating Male Luer	16.4 mL	50/Case
					
MRA102-84 ©	10	84 in. (213 cm)	Vented/Non-vented, Back Check Valve, Pre-pierced Y-site, (2) MicroClave® Clear Needle-free Y-sites	10.7 mL	50/Case
					
MRSA103-94	10	94 in. (239 cm)	Vented/Non-vented, Back Check Valve, Pre-pierced Y-site, (3) MicroClave Clear Needle-free Y-sites	17.5 mL	50/Case
					
MRSA10E-110	10	110 in. (279 cm)	Vented/Non-vented, Back Check Valve, (2) Pre-pierced Y-sites, (2) MicroClave Clear with Pre-attached Ext Set	21.1 mL	50/Case
					

## 15 Drops/mL

### General Primary Sets

REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
AA3101 ©	15	78 in. (198 cm)	Vented/Non-vented, Slide Clamp, Injection Y-site	15 mL	50/Case
					
158310 ©	15	83 in. (211 cm)	Vented/Non-vented, Pre-pierced Y-site	17 mL	50/Case
					

## 15 Drops/mL

### General Primary Sets

REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
158309 ©	15	83 in. (211 cm)	Vented/Non-vented, Sure-Lok® Needle-free Y-site	17 mL	50/Case
					
MRSA15-100	15	100 in. (254 cm)	Vented/Non-vented, MicroClave® Clear Needle-free Y-site	12.4 mL	50/Case
					
158306 ©	15	83 in. (211 cm)	Pre-pierced Y-site, Sure-Lok Needle-free Y-site	16 mL	50/Case
					
MRSA15-91	15	91 in. (231 cm)	Vented/Non-vented, Pre-pierced Y-site, MicroClave Clear Needle-free Y-site	11.4 mL	50/Case
					
159602 ©	15	96 in. (244 cm)	(2) Injection Y-sites	17.7 mL	50/Case
					
MRA15-100	15	100 in. (254 cm)	Vented/Non-vented, Pre-pierced Y-site, MicroClave Clear Needle-free Y-site	13 mL	50/Case
					
1510501 ©	15	105 in. (266 cm)	Vented/Non-vented, (2) Sure-Lok Needle-free Y-sites	20 mL	50/Case
					
158305 ©	15	89 in. (226 cm)	Pre-pierced Y-site, Sure-Lok Needle-free Y-site with Preattached Ext Set, Pinch Clamp	17.5 mL	50/Case
					
MRSA15E-88	15	88 in. (224 cm)	Vented/Non-vented, Pinch Clamp, Pre-pierced Y-site with Pre-attached Ext Set, Pinch Clamp, (2) MicroClave Clear Needle-free Y-sites	11.1 mL	50/Case
					
158307 ©	15	89 in. (226 cm)	Sure-Lok Needle-free Y-site, Pre-pierced Y-site with Preattached Ext Set, Sure-Lok Needle-free Y-site, Slide Clamp, Pre-pierced Y-site	18.7 mL	50/Case
					
MRSA15E-120	15	120 in. (305 cm)	Vented/Non-vented, Pre-pierced Y-site, MicroClave Clear Needle-free Y-site, Rotating Male Luer Lock, Drop-In Ext Set MicroClave Clear Y-site	15.8 mL	50/Case
					










# Gravity Administration IV Sets

## 15 Drops/mL

### Piggyback Sets with Back Check Valve

- Provides the flexibility to infuse more than one bag of medication or solution
- Featuring a back check valve to prevent fluid exchange or backflow into the primary infusion line

REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
MRSA15E-96	15	96 in. (244 cm)	Vented/Non-vented, Back Check Valve, (2) MicroClave® Clear Needle-free Y-sites, Pre-attached Ext Set, Slide Clamp	12.1 mL	50/Case
					
MRA15B-105	15	105 in. (267 cm)	Vented/Non-vented, Back Check Valve, (2) MicroClave Clear Needle-free Y-sites	13.2 mL	50/Case
					
AS3100 ©	15	110 in. (279 cm)	Vented/Non-vented, Back Check Valve, (2) Injection Y-sites	19.7 mL	50/Case
					
AS3101 ©	15	110 in. (279 cm)	Vented/Non-vented, Back Check Valve, Slide Clamp, (2) Sure-Lok™ Needle-free Y-sites, Roller Clamp	20.8 mL	50/Case
					
AS3102 ©	15	110 in. (279 cm)	Vented/Non-vented, Back Check Valve, Slide Clamp, (2) Pre-pierced Y-sites, Roller Clamp	20.4 mL	50/Case
					
AS3110A ©	15	110 in. (279 cm)	Vented/Non-vented, Back Check Valve, (3) Sure-Lok Needle-free Y-sites, Roller Clamp, Slide Clamp	22.3 mL	50/Case
					
AS3110B ©	15	110 in. (279 cm)	Vented/Non-vented, Back Check Valve, (3) Pre-pierced Y-sites, Slide Clamp	22 mL	50/Case
					



### 60 Drops/mL General Primary Sets

REF	DROPS/mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
607201 ©	60	72 in. (183 cm)	Injection Y-site	14 mL	50/Case
608304	60	89 in. (226 cm)	Non-vented, Roller Clamp (2), Sure-Lok® Needle-free Y-site with Pre-attached Ext Set with Pinch Clamp	17.5 mL	50/Case
608305 ©	60	89 in. (226 cm)	Pre-pierced Y-site, Sure-Lok Needle-free Y-site with Pre-attached Ext Set, Pinch Clamp	17.5 mL	50/Case
608306 ©	60	83 in. (211 cm)	Pre-pierced Y-site, Sure-Lok Needle-free Y-site	16 mL	50/Case
MRSA60-84	60	84 in. (213 cm)	Vented/Non-vented, Pre-pierced Y-site, (2) MicroClave® Clear Needle-free Y-sites	10.7 mL	50/Case
MRSA60E-88	60	88 in. (224 cm)	Vented/Non-vented, Pre-pierced Y-site with Pre-attached Ext Set, (2) MicroClave Clear Needle-free Y-sites	11.1 mL	50/Case
MRSA60-100	60	100 in. (254 cm)	Vented/Non-vented, MicroClave Clear Needle-free Y-site, Pre-pierced Y-site	12.5 mL	50/Case

### 60 Drops/mL Piggyback Sets with Back Check Valve

- Provides the flexibility to infuse more than one bag of medication or solution
- Featuring a back check valve to prevent fluid exchange or backflow into the primary infusion line

REF	DROPS/mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
MRA60B-92	60	92 in. (234 cm)	Vented/Non-vented, Back Check Valve, (2) MicroClave Clear Needle-free Y-sites, Slide Clamp	11.4 mL	50/Case
MRSA602-100	60	100 in. (254 cm)	Vented/Non-vented, Back Check Valve, (2) MicroClave Clear Needle-free Y-sites	12.5 mL	50/Case

### 60 Drops/mL Burette Sets

REF	DROPS/mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
MRS60B-106 ©	60	106 in. (269 cm)	150 mL Burette Set (w/MicroClave®, Shut-Off), (2) MicroClave Clear Needle-free Y-sites	12.7 mL	25/Case
MRS60B-117 ©	60	117 in. (297 cm)	24-hour 150 mL Burette Set (w/15 Micron Filter, Shut-off) with (3) MicroClave Clear Needle-free Y-sites	13.5 mL	25/Case

## AMSafe-3® Selectable Drip Chamber

### 10, 15 and 60 Drops/mL

AMSafe-3 is a 3-in-1 multi-drip chamber IV set that allows providers to select a 10, 15 or 60 drops/mL setting. Designed for efficiency where inventory reduction and cost containment are essential.

- Drop setting can be adjusted during medication delivery
- No need for clamps above and below chamber
- Includes male rotating luer lock
- Sterile, DEHP-free, not manufactured with natural rubber latex



REF	DROPS/mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
D38301 ©	10, 15 or 60	83 in. (210 cm)	Pre-pierced Y-site, Sure-Lok® Needle-free Y-site	27.7 mL	50/Case
D38302 ©	10, 15 or 60	91 in. (231 cm)	Pre-pierced Y-site, Sure-Lok Needle-free Y-site with Pre-attached Ext Set with Pre-pierced Y-site, Slide Clamp, Sure-Lok Needle-free Y-site	30.2 mL	50/Case
D38305 ©	10, 15 or 60	91 in. (231 cm)	Pre-pierced Y-site, Sure-Lok Needle-free Y-site, Back Check Valve with Pre-attached Ext Set with Pre-pierced Y-site, Slide Clamp, Sure-Lok Needle-free Y-site	30.5 mL	50/Case




## Secondary Piggyback Sets

REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
<b>AS3515</b> ©	15	40 in. (102 cm)	Vented/Non-vented, Roller Clamp, Hanger	10 mL	50/Case
<b>AS3500</b> ©	20	40 in. (102 cm)	Vented/Non-vented, Roller Clamp, Hanger	9.2 mL	50/Case




## Y-Type Blood Sets


REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
<b>MRSB-94</b> ©	10	94 in. (239 cm)	170 Micron Filter, Roller Clamps, Hand Pump, Slide Clamp, Roller Clamp, MicroClave® Y-site	57.2 mL	25/Case



<b>MRSB-102</b>	20	102 in. (259 cm)	200 Micron Filter, Roller Clamps, Hand Pump, MicroClave Y-site	20.2 mL	25/Case
-----------------	----	------------------	--	---------	---------



<b>MRSB-105</b>	10	105 in. (267 cm)	170 Micron Filter, Roller Clamps, Hand Pump, Pre-pierced Y-site, Slide Clamp, Roller Clamp, MicroClave Y-site, Pre-attached Ext Set	59.3 mL	25/Case
-----------------	----	------------------	---	---------	---------



## IV Extension Sets

Amsino's extension sets are available in various lengths and dimensions to accommodate clinical applications in healthcare settings. All extension sets are DEHP-free, not manufactured with natural rubber latex, sterile and include a rotating male luer lock unless otherwise noted.

### Standard Extension Sets

REF	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
<b>000606</b>	6 in. (15 cm)	Female Luer Lock, Pre-pierced Y-site, Sure-Lok® Needle-free Y-site, Slide Clamp	2.5 mL	50/Case



<b>MR4T07</b> ©	6 in. (15 cm)	MicroClave® Clear Needle-free Connector, T-connector	0.86 mL	50/Case
-----------------	---------------	--	---------	---------



<b>AE3107</b> ©	7 in. (18 cm)	Removeable Sure-Lok Needle-free Connector, Slide Clamp	1.76 mL	50/Case
-----------------	---------------	--	---------	---------



<b>MR4107</b> ©	7 in. (18 cm)	Removable MicroClave Clear Needle-free Connector, Slide Clamp	0.69 mL	50/Case
-----------------	---------------	---	---------	---------



<b>MR4207</b> ©	7 in. (18 cm)	Bonded MicroClave Clear Needle-free Connector, Slide Clamp	0.69 mL	50/Case
-----------------	---------------	--	---------	---------



<b>AE3108</b> ©	8 in. (20 cm)	Removeable Sure-Lok Needle-free Connector, Pinch Clamp	1.76 mL	100/Case
-----------------	---------------	--	---------	----------



<b>MR4R10</b> ©	10 in. (25 cm)	Removable MicroClave Clear Needle-free Connector, Pinch Clamp	0.98 mL	50/Case
-----------------	----------------	---	---------	---------



<b>003201</b> ©	27 in. (69 cm)	Female Luer Lock, (2) Slide Clamps, (2) Pre-pierced Y-sites	4.9 mL	50/Case
-----------------	----------------	---	--------	---------



<b>003202</b>	27 in. (69 cm)	Female Luer Lock, (2) Slide Clamps, (2) Sure-Lok Needle-free Y-sites	5.3 mL	50/Case
---------------	----------------	--	--------	---------



<b>MR4Y32</b> ©	32 in. (81 cm)	(2) MicroClave Clear Needle-free Connectors, (2) Slide Clamps	5.02 mL	50/Case
-----------------	----------------	---	---------	---------



## IV Extension Sets

### Microbore Extension Sets

REF	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
<b>AE0106</b> ©	6 in. (15 cm)	Removeable Sure-Lok® Needle-free Connector, Slide Clamp	0.63 mL	50/Case
<b>AE1106</b> ©	6 in. (15 cm)	Female Luer Lock, Slide Clamp, Fixed Male Luer Lock	0.5 mL	50/Case
<b>AE0107</b> ©	7 in. (18 cm)	Female Luer Lock, Slide Clamp, Sure-Lok Needle-free Y-site	0.9 mL	50/Case
<b>AE1112</b> ©	12 in. (31 cm)	Female Luer Lock, Slide Clamp	0.7 mL	50/Case
<b>AE1138</b> ©	38 in. (97 cm)	Female Luer Lock, Slide Clamp	1.0 mL	50/Case

### Smallbore Extension Sets

REF	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
<b>MR4307</b> ©	7 in. (18 cm)	Bonded MicroClave® Clear Needle-free Connector, Slide Clamp	0.24 mL	50/Case
<b>MR4B07</b> ©	7 in. (18 cm)	Bifuse Ext Set with (2) MicroClave Clear Needle-free Connector, (2) Slide Clamps	0.35 mL	50/Case
<b>MR4TF04</b> ©	4 in. (10 cm)	Trifuse Ext Set with (3) MicroClave Clear Needle-free Connectors, (3) Slide Clamps	0.4 mL	50/Case

### Flow Controller Extension Sets

Amsino's IV extension sets with flow control devices are designed to control IV fluid and medications using a 'regulator dial' to control the rate of infusion. Includes a rotating male luer lock.

REF	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
<b>AFS107</b> ©	20 in. (51 cm)	Flow Controller, Injection Y-site, Slide Clamp	3.5 mL	50/Case
<b>MRSFC-19</b> ©	19 in. (48 cm)	Flow Controller, MicroClave Clear Needle-free Connector, Clamp	3.2 mL	50/Case

## Pressure Infusion Extension Sets

Amsino's single-patient-use pressure infusion IV extension set is made to withstand up to 400 psig for rapid administration of medication or fluids in an emergency. Includes a rotating male luer lock.

REF	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
<b>MR4P07</b> ©	7 in. (18 cm)	Pressure Infusion (400psig) Ext Set with MicroClave® Clear Needle-free Connector, Purple Slide Clamp	0.43 mL	50/Case

## Needle-Free Connectors

Prevent needle stick injuries and reduce bacterial contamination with needle-free connectors.

- Swabbable
- MRI-safe
- Not limited to 24-hour use
- Sterile, DEHP-free, not manufactured with natural rubber latex

REF	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
<b>MR4001</b>	MicroClave Clear Needle-free Pressure Infusion (400psig) Connector	0.04 mL	100/Case
<b>AY0200</b> ©	AMSafe® Needle-free Connector	0.04 mL	100/Case

## Bag Spike

REF	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
<b>MBSS001</b>	MicroClave Needle-free Connector IV Bag Spike	0.34 mL	50/Case

## Needle-Free Vial Adapter

REF	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
<b>MB4V01</b>	MicroClave Clear Needle-free Connector Single Vial Adapter	0.16 mL	50/Case



# IV Pump Administration Sets




## Sure-Lok® Needle-Free IV Pump Administration Sets

- FDA cleared IV pump sets for use with Baxter Spectrum and Flo-Gard® 6201 and 6301 series pumps
- Sure-Lok Needle-free connector and Pre-Pierced Y-injection site options
- 0.22 micron filter option available
- TUBETECH® technology provides +/- 5% accurate fluid delivery over 24 hour timeframe
- DEHP-free, non-pyrogenic, not manufactured with natural rubber latex

REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
FP092	10	92 in. (234 cm)	Roller Clamp, Rotating Male Luer Lock with Retractable Collar	14 mL	50/Case
					
FP092Y1	10	92 in. (234 cm)	1 Needle-free Valve, Roller Clamp, 1 Slide Clamp, Rotating Male Luer Lock with Retractable Collar	14 mL	50/Case
					
FP100BY2	10	100 in. (254 cm)	Back check Valve, 2 Needle-free Valves, Roller Clamp, 2 Slide Clamps, Rotating Male Luer Lock with Retractable Collar	16 mL	50/Case
					
FP109BY3	10	109 in. (277 cm)	Back check Valve, 3 Needle-free Valves, Roller Clamp, 2 Slide Clamps, Rotating Male Luer Lock with Retractable Collar	18 mL	50/Case
					
FP169BY3	60	109 in. (277 cm)	Needle-free Valve, Roller Clamp, 2 Slide Clamps, Back check Valve, Rotating Male Luer Lock with Retractable Collar	18 mL	50/Case
					
FP106FY2	10	106 in. (269 cm)	Back check Valve, 2 Needle-free Valves, Roller Clamp, 2 Slide Clamps, 0.22 Micron Filter, Rotating Male Luer lock with Retractable Collar	22 mL	50/Case
					
FP100FY1	10	100 in. (254 cm)	1 Needle-free Valve, Roller Clamp, 2 Slide Clamps, 0.22 Micron Filter, Rotating Male Luer lock with Retractable Collar	22 mL	50/Case
					

## Pre-Pierced Y-Injection IV Pump Administration Sets

- 510(k) cleared IV pump sets for use with Baxter Spectrum and Flo-Gard® 6201 and 6301 series pumps

REF	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOLUME	QTY/UNIT
FP100BI2	10	100 in. (254 cm)	Back check Valve, 2 Pre-pierced Injection Y-sites, Roller Clamp, 2 Slide Clamps, Rotating Male Luer Lock with Retractable Collar	16 mL	50/Case
					
FP109BI3	10	109 in. (277 cm)	Back check Valve, 3 Pre-pierced Injection Y-Sites, Roller Clamp, 2 Slide Clamps, Rotating Male Luer Lock with Retractable Collar	18 mL	50/Case
					
FP106FI2	10	106 in. (269 cm)	Back check Valve, 2 Pre-pierced Injection Y-sites, Roller Clamp, 2 Slide Clamps, 0.22 Micron Filter, Rotating Male Luer lock with Retractable Collar	22 mL	50/Case
					

## AMSafe® Prefilled Normal Saline Flush Syringe

- 0.9% Sodium Chloride Injection, USP
- 3 mL, 5 mL, 10 mL and 20 mL
- Sterile Field Ready (SFR) and Sterile Fluid Path packaging configurations



For more details and ordering information, please refer to page 55.

### Solutions for Infection Prevention

AMSure Urological Care product portfolio provides solutions to help improve patient safety and infection prevention. Our product selection includes:

- Foley Catheter Trays & Collection ..... 22
  - Pre-Connected Foley Catheter Urine Drainage Bag Trays..... 22
  - Add-a-Foley Catheter Trays ..... 23
  - Pre-Connected Foley Catheter Urine Meter Kits ..... 24
  - Add-a-Foley Urine Meter Trays ..... 25
- Foley Insertion Trays..... 25
- Urinary Drainage Bags ..... 26
- Urine Meter ..... 27
  - ACCU-MEASURE® Pediatric Urine Meter ..... 27
- Foley Catheters ..... 28
  - Hydrophilic Latex Foley Catheters ..... 28
  - Silicone-Coated Latex Foley Catheters ..... 29
  - 100% Silicone Foley Catheters ..... 30
- Intermittent Urinary Catheters ..... 31
  - Urethral Insertion Trays and Kits ..... 31
  - Closed System Urethral Drainage Kits ..... 32
  - Self-Catheterization Kits with R-Polished™ Eye Catheters ..... 33
  - PVC Intermittent Urethral Catheter with R-Polished™ Eyes..... 34
  - PVC Intermittent Urethral Catheter ..... 35
  - Red Rubber Latex Urethral Catheter ..... 35
- Irrigation Sets..... 36
  - Cystoscopy & TUR Irrigation Sets ..... 36
- Ancillary & Patient Care Products..... 37
  - Urinary Leg Bags ..... 37
  - Extension Tube for Leg Bags ..... 37
  - Pediatric Urine Collector ..... 37
  - Irrigation Trays and Kits ..... 38
  - Syringes ..... 39
  - Specimen Containers ..... 40

### Foley Catheter Trays & Collection



AS89316

#### AMSure® Pre-Connected Foley Catheter Urine Drainage Bag Trays

AMSure® Pre-Connected Drainage Bag Foley Trays improve simplicity and efficiency of care. The catheter is pre-attached to a low-profile drainage bag which helps to minimize potential contamination. (HCPCS codes: See product chart below)

- Tray includes all necessary components for aseptic catheter insertion
- Foley catheter connected with tamper evident seal to prevent disconnection and maintain sterile fluid path
- Tray improves clinical efficacy through cost-effective inventory management
- Gloves and povidone iodine solution help facilitate aseptic technique and reduce the risk of infection
- Needle-free urinary sampling port reduces the risk of needle stick injuries
- Trays are available with 100% silicone, silicone-coated latex, or hydrophilic-coated latex Foley catheters to suit patient needs
- Sterile

**Trays include:** 2000 mL low-profile urinary drainage bag, pre-connected 16 in. Foley catheter with 5 cc balloon, waterproof and fenestrated drapes, Nitrile powder-free gloves, (5) Rayon balls, 10% Povidone Iodine solution, forceps, lubricating jelly, syringe prefilled with 10 cc sterile water, 2 oz. specimen container with a patient label, central supply room wrap and Tyvek header bag

CATHETER TYPE	REF	FOLEY SIZE	BALLOON SIZE	QTY/UNIT
100% Silicone (HCPCS Code A4315)	AS89314S	14 Fr.	5 cc	10/Case
	AS89316S	16 Fr.	5 cc	10/Case
	AS89318S	18 Fr.	5 cc	10/Case
Silicone Coated Latex (HCPCS Code A4314)	AS89314	14 Fr.	5 cc	10/Case
	AS89316	16 Fr.	5 cc	10/Case
	AS89318	18 Fr.	5 cc	10/Case
Hydrophilic Coated Latex (HCPCS Code A4314)	AS89314H	14 Fr.	5 cc	10/Case
	AS89316H	16 Fr.	5 cc	10/Case
	AS89318H	18 Fr.	5 cc	10/Case



## Foley Catheter Trays & Collection



AS89300

### AMSure® Add-a-Foley Catheter Trays

AMSure Add-a-Foley Catheter Trays improve simplicity and efficiency of care by including the urine collection device and necessary supplies to safely insert a Foley catheter. (HCPCS code A4354)

- Needle-free sampling port on urinary bag reduces the risk of needle stick injuries
- Gloves and povidone iodine swab sticks help maintain sterility and minimized risk of infection
- Sterile, not manufactured with natural rubber latex

Trays include: 2000 mL low-profile urinary drainage bag, waterproof and fenestrated drapes, Nitrile powder-free gloves, (5) Rayon balls, 10% povidone iodine solution, forceps, lubricating jelly, 2 oz. specimen container with patient label, central supply room wrap and Tyvek header bag

REF	DESCRIPTION	PRE-FILLED SYRINGE	QTY/UNIT
<b>AS89300</b>	Add-a-Foley Tray 2000 mL Urine Drainage Bag	10 cc	10/Case
<b>AS89330</b> ©	Add-a-Foley Tray 2000 mL Urine Drainage Bag	30 cc	10/Case

## AMSure® Pre-Connected Foley Catheter Urine Meter Kits

The AMSure pre-connected urine meter trays offer the benefit of a closed system along with the accuracy of a meter to measure exact urine output. These truly closed systems provide all the necessary supplies for an aseptic catheter insertion to promote the US Centers for Disease Control (CDC) recommendations to help prevent catheter associated urinary tract infections (CAUTIs). (HCPCS codes: See product chart below)

Trays include: 400 mL Urine Meter Preconnected with 5 cc Foley Catheter, waterproof and fenestrated drapes, (2) Nitrile gloves, 10% Povidone Iodine solution, (5) Rayon balls, forceps, lubricating jelly, syringe pre-filled with 10 cc sterile water, 2 oz. specimen container with a patient label, central supply room wrap, Tyvek header bag



AS898145

CATHETER TYPE	REF	FOLEY SIZE	BALLOON SIZE	QTY/UNIT
100% Silicone (HCPCS Code A4315)	<b>AS89814S</b> ©	14 Fr.	5 cc	10/Case
	<b>AS89816S</b> ©	16 Fr.	5 cc	10/Case
	<b>AS89818S</b> ©	18 Fr.	5 cc	10/Case
Silicone Coated Latex (HCPCS Code A4314)	<b>AS89814</b> ©	14 Fr.	5 cc	10/Case
	<b>AS89816</b> ©	16 Fr.	5 cc	10/Case
	<b>AS89818</b> ©	18 Fr.	5 cc	10/Case
Hydrophilic Coated Latex (HCPCS Code A4314)	<b>AS89814H</b> ©	14 Fr.	5 cc	10/Case
	<b>AS89816H</b> ©	16 Fr.	5 cc	10/Case
	<b>AS89818H</b> ©	18 Fr.	5 cc	10/Case

## Foley Catheter Trays & Collection

### AMSure® Add-a-Foley Urine Meter Trays

The AMSure Add-a-Foley urine meter kits provide all the supplies necessary for an aseptic insertion. They also offer the caregiver the flexibility to choose the best Foley catheter for the patient and the situation, along with providing a urine meter to obtain the most accurate urine output measurements. (HCPCS Code A4354)



AS89800

- Sterile, not manufactured with natural rubber latex

Trays include: 400 ml Urine Meter Preconnected with 5 cc Foley Catheter, waterproof and fenestrated drapes, (2) Nitrile gloves, 10% Povidone Iodine solution, (5) Rayon balls, forceps, lubricating jelly, syringe prefilled with 10 cc sterile water, 2 oz. specimen container with a patient label, central supply room wrap, Tyvek header bag

REF	DESCRIPTION	PRE-FILLED SYRINGE	QTY/UNIT
AS89800	Add-a-Foley Tray 400 mL Urine Meter with 2600 mL Urine Drainage Bag	10 cc	10/Case
AS89830	Add-a-Foley Tray 400 mL Urine Meter with 2600 mL Urine Drainage Bag	30 cc	10/Case

### AMSure® Foley Insertion Trays

AMSure Foley Insertion Trays improve simplicity and efficiency of care by providing the basic supplies required to safely insert a Foley catheter. (HCPCS code A4310)



AS880

- Deep tray basin designed to prevent contamination
- Gloves and povidone iodine swab sticks help maintain sterility and minimized risk of infection
- Sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex

Trays include: 1000 mL collection basin, (3) Povidone-iodine swab sticks or benzalkonium chloride prep pad, vinyl powder-free gloves, waterproof underpad 13 in. x 18 in. (33.02 cm x 45.72 cm), lubrication jelly, fenestrated drape

SWAB STICK	REF	PRE-FILLED SYRINGE	QTY/UNIT
PVP Swab Sticks	AS880 ©	10 cc	10/Case
	AS890 ©	30 cc	10/Case
BZK Swab Sticks	AS880K ©	10 cc	10/Case
	AS890K ©	30 cc	10/Case

## AMSure® Urinary Drainage Bags

AMSure urine drainage bags are available in a variety of different configurations to meet the clinicians' needs to best provide patient safety and comfort while helping to minimize infections. All drainage bags include a sheet clamp for the positioning and securing of the drainage tube. The drainage tube has an ID of 1 1/32" (0.87 cm) for unrestricted flow and are 4 feet (1.2 M) in length, terminating with a standard terraced catheter fitting for a secure and leak-proof connection. Each bag includes a storage pocket for the emptying tube. (HCPCS code A4358)

- 2,000 cc capacity graduated in 200 cc increments
- Anti-reflux chambers help improve patient safety by restricting urine backflow
- Leak-proof seams ensure fluid containment
- Needle-free sampling port reduces the risk of needle stick injuries
- Air vents help to maintain consistent urine flow, delivering safer patient outcomes
- Sterile fluid pathway, BPA-free, not manufactured with natural rubber latex



AS322 Teardrop



AS326 Oval



AS32600 Low-Profile

REF	BAG SHAPE	ANTI-REFLUX TYPE	NEEDLE-FREE SAMPLING PORT	DRAIN TYPE	HANGER TYPE	QTY/UNIT
AS312 ©	Diamond	Drip Device	Pre-Pierced Blunt Cannula	T-Tap	Double-Hook & Rope	20/Case
AS322 ©	Teardrop	Drip Device	Pre-Pierced Blunt Cannula	T-Tap	Double-Hook & Rope	20/Case
AS332	Teardrop	Small Chamber	Pre-Pierced Blunt Cannula	T-Tap	Double-Hook & Rope	20/Case
AS302 ©	Pear	Flutter Valve	Pre-Pierced Blunt Cannula	T-Tap	Single-Hook	20/Case
AS32200 ©	Low-Profile	Large Chamber	Pre-Pierced Blunt Cannula	T-Tap	Swivel Single-Hook & Rope	20/Case
AS32600* ©	Low-Profile	Large Chamber	Male Luer-Lock or Luer-Slip	T-Tap	Swivel Single-Hook & Rope	20/Case
AS326* ©	Oval	Drip Device	Pre-Pierced Blunt Cannula	Metal Clamp with Doghouse Storage	Double-Hook & Rope	20/Case

\*100% sterile packed in peel pouch



## Foley Catheter Trays & Collection

### AMSure® Urine Meter

AMSure Urine Meter Drain Bags are designed to provide all caregivers the ability to easily and accurately measure urine output. The meter device is pre-attached to a urine drainage bag, for convenient draining of the meter while still obtaining incremental urine measurements. (HCPCS code A4358)

- **Anti-Reflux Chamber helps improve patient safety by reducing urine backflow**
- **Twist valve drains at the bottom of the urine meter for easy drainage**
- **Leak-proof seams ensure fluid containment**
- **Needle-free sampling port reduces the risk of needle stick injuries**
- **Air vents help to maintain consistent urine flow delivering safer patient outcomes**
- **Sterile, not manufactured with natural rubber latex**

REF	DESCRIPTION	QTY/UNIT
<b>AS328</b> ©	400 mL Urine Meter with 2600 mL Urine Drainage Bag	10/Case



### ACCU-MEASURE® Pediatric Urine Meter

2000 mL Bag with 150 mL Chamber

ACCU-MEASURE Pediatric Urine Meter is specifically designed for pediatric patients. With 1 mL increments, the 150 mL meter chamber's white background provides accurate and precise output measurements. The needle-free sample port reduces the risk of needle stick injuries, helping keep nurses safe while making their job easier. (HCPCS code A4358)

- **Double outlet valves – Twist valve for taking sample, T-valve for emptying the bag**
- **Needle-free sampling port with luer lock access**
- **Leak-proof seams ensure fluid containment**
- **Meter bottle capacity 150 mL**
- **Adjustable bed sheet clamp for securing drainage bag tubing**
- **Sterile, BPA-free, not manufactured with natural rubber latex**

REF	DESCRIPTION	QTY/UNIT
<b>AS327100</b> ©	ACCU-MEASURE® Pediatric Urine Meter 2000 mL Bag with 150 mL Chamber	30/Case



## Foley Catheters

### AMSure® Hydrophilic Latex Foley Catheters

(HCPCS code A4338)

- **Hydrophilic coating provides for lubricious surface when contacting moisture to help reduce irritation of the urethra**
- **Smooth eyes for improved patient comfort**
- **Symmetrical balloon provides a tight seal at the neck of the urethra**
- **Sterile, DEHP-free, BPA-free**



BALLOON SIZE	REF	FOLEY SIZE	QTY/UNIT
5 cc	<b>AS41014H</b> ©	14 Fr. x 16 in. (40.6 cm)	10/Box
	<b>AS41016H</b> ©	16 Fr. x 16 in. (40.6 cm)	10/Box
	<b>AS41018H</b> ©	18 Fr. x 16 in. (40.6 cm)	10/Box
30 cc	<b>AS42014H</b>	14 Fr. x 16 in. (40.6 cm)	10/Box
	<b>AS42016H</b> ©	16 Fr. x 16 in. (40.6 cm)	10/Box
	<b>AS42018H</b> ©	18 Fr. x 16 in. (40.6 cm)	10/Box

# Foley Catheters

## AMSure® Silicone-Coated Latex Foley Catheters

(HCPCS code A4338)


- Silicone coating helps reduce urethral irritation
- Smooth beveled symmetrical eyes minimize tissue damage
- Large oval eyes aid in maximum drainage
- Firm catheter tip designed to improve ease of insertion
- Symmetrical balloon provides a tight seal at the neck of the urethra
- Individually packaged, sterile, DEHP-free, BPA-free



## AMSure® 100% Silicone Foley Catheters

(HCPCS code A4344)

- Eliminates latex allergy concerns
- Smooth eyes improve patient comfort
- Reinforced tip provides a stylet for optimized insertion efficacy
- 100% silicone material has non-adherent characteristics and reduces encrustation and calcifications
- Symmetrical balloon provides a tight seal at the neck of the urethra
- Individually packaged, sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex



BALLOON SIZE	REF	FOLEY SIZE	QTY/UNIT
3 cc	AS40008 ©	8 Fr. x 10.6 in. (27cm)	10/Box
	AS40010 ©	10 Fr. x 10.6 in. (27cm)	10/Box
5 cc	AS41012 ©	12 Fr. x 16 in. (40.6cm)	10/Box
	AS41014 ©	14 Fr. x 16 in. (40.6cm)	10/Box
	AS41016 ©	16 Fr. x 16 in. (40.6cm)	10/Box
	AS41018 ©	18 Fr. x 16 in. (40.6cm)	10/Box
	AS41020 ©	20 Fr. x 16 in. (40.6cm)	10/Box
	AS41022 ©	22 Fr. x 16 in. (40.6cm)	10/Box
	AS41024 ©	24 Fr. x 16 in. (40.6cm)	10/Box
	AS41026 ©	26 Fr. x 16 in. (40.6cm)	10/Box
	AS41028 ©	28 Fr. x 16 in. (40.6cm)	10/Box
	AS41030 ©	30 Fr. x 16 in. (40.6cm)	10/Box
30 cc	AS42012 ©	12 Fr. x 16 in. (40.6cm)	10/Box
	AS42014 ©	14 Fr. x 16 in. (40.6cm)	10/Box
	AS42016 ©	16 Fr. x 16 in. (40.6cm)	10/Box
	AS42018 ©	18 Fr. x 16 in. (40.6cm)	10/Box
	AS42020 ©	20 Fr. x 16 in. (40.6cm)	10/Box
	AS42022 ©	22 Fr. x 16 in. (40.6cm)	10/Box
	AS42024 ©	24 Fr. x 16 in. (40.6cm)	10/Box
	AS42026 ©	26 Fr. x 16 in. (40.6cm)	10/Box
	AS42028 ©	28 Fr. x 16 in. (40.6cm)	10/Box
AS42030 ©	30 Fr. x 16 in. (40.6cm)	10/Box	

BALLOON SIZE	REF	FOLEY SIZE	QTY/UNIT
5 cc	AS41012S ©	12 Fr. x 16 in. (40.6cm)	10/Box
	AS41014S ©	14 Fr. x 16 in. (40.6cm)	10/Box
	AS41016S ©	16 Fr. x 16 in. (40.6cm)	10/Box
	AS41018S ©	18 Fr. x 16 in. (40.6cm)	10/Box
	AS41020S ©	20 Fr. x 16 in. (40.6cm)	10/Box
	AS41022S ©	22 Fr. x 16 in. (40.6cm)	10/Box
	AS41024S ©	24 Fr. x 16 in. (40.6cm)	10/Box
	AS41026S ©	26 Fr. x 16 in. (40.6cm)	10/Box
	AS41028S ©	28 Fr. x 16 in. (40.6cm)	10/Box
30 cc	AS42012S ©	12 Fr. x 16 in. (40.6cm)	10/Box
	AS42014S ©	14 Fr. x 16 in. (40.6cm)	10/Box
	AS42016S ©	16 Fr. x 16 in. (40.6cm)	10/Box
	AS42018S ©	18 Fr. x 16 in. (40.6cm)	10/Box
	AS42020S ©	20 Fr. x 16 in. (40.6cm)	10/Box



# Intermittent Urinary Catheters



AS870

## AMSure® Urethral Insertion Trays and Kits

AMSure® Urethral Trays improve simplicity and efficiency of care by including all components needed to perform short-term aseptic catheterization. (HCPCS code A4353)

- **Gloves and povidone iodine swab sticks help maintain aseptic technique and reduce the risk of infection**

### PVC Catheter Tray

- **Staggered eyes help maintain tip strength during insertion and promote efficient drainage**
- **Short-term catheter use helps decrease infections**
- **Sterile, BPA-free, not manufactured with natural rubber latex**

### Red Rubber Catheter Trays

- **Smooth tip and quality rubber latex promote patient comfort during insertion**
- **Clearly marked French size on funnel promotes ease of selection**
- **Sterile, DEHP-free, BPA-free**

Trays include: 1000 mL outer basin tray, 16 in. urethral catheter, waterproof and fenestrated drapes, vinyl powder-free gloves, (3) 10% povidone iodine swab sticks, lubricating jelly, 4 oz. specimen container

REF	URETHRAL CATHETER TYPE	QTY/UNIT
<b>AS870</b> ©	14 Fr PVC*	20/Case
<b>AS87114</b> ©	14 Fr Red Rubber Latex	20/Case
<b>AS87118</b> ©	18 Fr Red Rubber Latex	20/Case

\*Not manufactured with natural rubber latex.

## AMSure® Closed System Urethral Drainage Kits

AMSure® Pre-connected Drainage Bag Urethral Catheter Kits improve simplicity and efficiency of care. The catheter is pre-attached to the drainage bag which helps to minimize contamination. (HCPCS code A4353)

- **Gloves and povidone iodine swab sticks help maintain sterility and reduce the risk of infection**
- **Sterile**

### Red Rubber Catheter Kits

- **Smooth tip and quality rubber latex promote patient comfort during insertion**
- **Clearly marked French size on funnel promotes ease of selection**
- **Sterile, DEHP-free, BPA-free**

### PVC Catheter Kits

- **Staggered eyes help maintain tip strength during insertion and promote efficient drainage**
- **Short-term catheter use helps decrease infections**
- **Sterile, BPA-free, not manufactured with natural rubber latex**



AS89214R

Kit includes: 1200 mL urinary drainage bag, 16 in. (40.64 cm) urethral catheter, CSR wrap, waterproof and fenestrated drapes, (2) Nitrile powder-free gloves, (3) Povidone-iodine swab sticks, lubricating jelly, 2 oz. specimen container with ID label

SWAB STICK	REF	FOLEY SIZE	FOLEY MATERIAL	QTY/UNIT
PVP Swab Sticks	<b>AS89214R</b> ©	14 Fr. x 16 in. (40.6cm)	Red Rubber (Contains Latex)	20/Case
	<b>AS89214</b> ©	14 Fr. x 16 in. (40.6cm)	PVC Latex- Free	20/Case
BZK Swab Sticks	<b>AS89214RK</b> ©	14 Fr. x 16 in. (40.6cm)	Red Rubber (Contains Latex)	20/Case
	<b>AS89216RK</b> ©	16 Fr. x 16 in. (40.6cm)	Red Rubber (Contains Latex)	20/Case

# Intermittent Urinary Catheters



AS85010

## AMSure® Self-Catheterization Kits with R-POLISHED™ Eye Catheters

AMSure Self-Catheterization kits with R-Polish eyelets technology allow the patient, even with limited dexterity, to perform incontinence care with confidence. R-Polish technology uses radio frequency to create ultra-smooth eyelets for optimal patient comfort.

The pre-lubricated, closed system with introducer tip helps to reduce potential for infections. The kits contain all the necessary components needed to perform aseptic self-catheterization. (HCPCS code A4353)

- Advanced R-Polished eyelets enhance patient comfort during catheter insertion and removal
- Introducer tip reduces the potential for urinary tract infection (UTI)
- Collection bag with built-in thumb holes allows for ease of use
- Sterile, BPA-free, not manufactured with natural rubber latex

Kit includes: Pre-lubricated 16 in. (40.64 cm) long PVC urethral catheter, 1100 mL collection bag with insertion tip and protector cap, (3) Povidone-iodine swab sticks and benzalkonium chloride prep pad, vinyl powder-free gloves, waterproof underpad 13 in. x 18 in. (33.02 cm x 45.72 cm), waste disposal bag, gauze pad 4 in. x 4 in. (10.16 cm x 10.16 cm)

SWAB STICK	REF	CATHETER SIZE	QTY/UNIT
PVP Swab Sticks	AS85008 ©	8 Fr. x 16 in. (40.6cm)	50/Case
	AS85010 ©	10 Fr. x 16 in. (40.6cm)	50/Case
	AS85012 ©	12 Fr. x 16 in. (40.6cm)	50/Case
	AS85014 ©	14 Fr. x 16 in. (40.6cm)	50/Case
	AS85016 ©	16 Fr. x 16 in. (40.6cm)	50/Case
BZK Swab Sticks	AS85008K ©	8 Fr. x 16 in. (40.6cm)	50/Case
	AS85010K ©	10 Fr. x 16 in. (40.6cm)	50/Case
	AS85012K ©	12 Fr. x 16 in. (40.6cm)	50/Case
	AS85014K ©	14 Fr. x 16 in. (40.6cm)	50/Case
	AS85016K ©	16 Fr. x 16 in. (40.6cm)	50/Case

## AMSure® PVC Intermittent Urethral Catheter with R-POLISHED™ Eyes

AMSure PVC Intermittent Catheters with R-Polish technology are designed for short-term use. R-Polish technology uses radio frequency to create ultra-smooth eyelets for optimal patient comfort. The intermittent catheters are available in various French sizes as well as female, pediatric and male lengths. (HCPCS code A4351)

- Advanced R-Polished eyelets enhance patient comfort during catheter insertion and removal
- Staggered eyes help maintain tip strength during insertion and promote efficient drainage
- Designed to minimize trauma and reduce risk of infection
- Sterile, BPA-free, not manufactured with natural rubber latex



FEATURES	REF	DESCRIPTION	QTY/UNIT
Female	AS960610	10 Fr, 6 in. (15.24 cm)	50/Case
	AS960612	12 Fr, 6 in. (15.24 cm)	50/Case
	AS960614	14 Fr, 6 in. (15.24 cm)	50/Case
	AS960616	16 Fr, 6 in. (15.24 cm)	50/Case
Pediatric	AS961010	10 Fr, 10 in. (25.40 cm)	50/Case
Male	AS961610	10 Fr, 16 in. (40.64 cm)	50/Case
	AS961612	12 Fr, 16 in. (40.64 cm)	50/Case
	AS961614	14 Fr, 16 in. (40.64 cm)	50/Case
	AS961616	16 Fr, 16 in. (40.64 cm)	50/Case
Coudé Tip	AS861610C	10 Fr, 16" (40.65 cm)	30/Case
	AS861612C	12 Fr, 16" (40.65 cm)	30/Case
	AS861614C	14 Fr, 16" (40.65 cm)	30/Case
	AS861618C	18 Fr, 16" (40.65 cm)	30/Case



## Intermittent Urinary Catheters

### AMSure® PVC Intermittent Urethral Catheter

AMSure PVC Intermittent Urethral Catheters provide a smooth surface firm enough for insertion yet soft enough to be slightly pliable. The urethral catheters are designed for intermittent use and available in various French sizes as well as female, pediatric and male lengths. (HCPCS code A4351)

- Staggered eyes help maintain tip strength during insertion and promote efficient drainage
- Short-term catheter use helps decrease infections
- Clear shaft provides easy viewing as you measure
- PVC material eliminates latex allergy concerns
- Individually packaged, sterile, BPA-free, not manufactured with natural rubber latex



FEATURES	REF	DESCRIPTION	QTY/UNIT
Female	AS860608 ©	8 Fr, 6 in. (15.24 cm)	50/Case
	AS860610 ©	10 Fr, 6 in. (15.24 cm)	50/Case
	AS860612 ©	12 Fr, 6 in. (15.24 cm)	50/Case
	AS860614 ©	14 Fr, 6 in. (15.24 cm)	50/Case
	AS860616 ©	16 Fr, 6 in. (15.24 cm)	50/Case
Pediatric	AS861006 ©	6 Fr, 10 in. (25.40 cm)	50/Case
	AS861008 ©	8 Fr, 10 in. (25.40 cm)	50/Case
	AS861010 ©	10 Fr, 10 in. (25.40 cm)	50/Case
Male	AS861608 ©	8 Fr, 16 in. (40.64 cm)	50/Case
	AS861610 ©	10 Fr, 16 in. (40.64 cm)	50/Case
	AS861612 ©	12 Fr, 16 in. (40.64 cm)	50/Case
	AS861614 ©	14 Fr, 16 in. (40.64 cm)	50/Case
	AS861616 ©	16 Fr, 16 in. (40.64 cm)	50/Case
	AS861618 ©	18 Fr, 16 in. (40.64 cm)	50/Case

### AMSure Red Rubber Urethral Catheter

AMSure Red Rubber Urethral Catheters are pliable and designed for intermittent use. The urethral catheters are available in various French sizes as well as female and male lengths. (HCPCS code A4351)

- Smooth tip and quality rubber latex promote patient comfort during insertion
- Clearly marked French size on funnel promotes ease of selection
- Individually packaged, sterile, DEHP-free




REF	DESCRIPTION	QTY/UNIT
AS44008 ©	8 Fr, 11 in. (27.94 cm)	100/Case
AS44010 ©	10 Fr, 11 in. (27.94 cm)	100/Case
AS44012 ©	12 Fr, 16 in. (40.64 cm)	100/Case
AS44014 ©	14 Fr, 16 in. (40.64 cm)	100/Case
AS44016 ©	16 Fr, 16 in. (40.64 cm)	100/Case
AS44018 ©	18 Fr, 16 in. (40.64 cm)	100/Case
AS44020 ©	20 Fr, 16 in. (40.64 cm)	100/Case
AS44022 ©	22 Fr, 16 in. (40.64 cm)	100/Case
AS44024 ©	24 Fr, 16 in. (40.64 cm)	100/Case
AS44028 ©	28 Fr, 16 in. (40.64 cm)	100/Case

## Irrigation Sets

### AMSure® Cystoscopy & TUR Irrigation Sets

The AMSure Cystoscopy and TUR Irrigation sets are used for continuous bladder irrigation after prostate or bladder surgeries. These are available with one, two, or four leads to meet a variety of needs. (HCPCS code A4351)

- Robert Clamp(s)
- Drip Chamber(s)
- Sterile, not manufactured with natural rubber latex



REF	TYPE	QTY/UNIT
ABL400 ©	Single Lead	20/Case

REF	TYPE	QTY/UNIT
ABL402 ©	Two Leads	20/Case



REF	TYPE	QTY/UNIT
ABL404 ©	Four Leads	20/Case

## Ancillary & Patient Care Products

### Urinary Leg Bags

(HCPCS code A4358)

- Adjustable, soft, breathable pre-attached elastic leg straps for patient comfort
- Anti-reflux flutter valve improves patient safety by reducing urine backflow
- Leak-proof seams ensure fluid containment
- Sterile fluid pathway
- BPA-free, not manufactured with natural rubber latex

SIZE	REF	DRAIN PORT STYLE	18 in. (45.72 cm) EXTENSION TUBE WITH ADAPTER & PROTECTIVE CAP	QTY/UNIT
600 mL	AS306N ©	Push-Pull		48/Case
	AS308N ©	Twist-Turn		48/Case
	LBE306N	Push-Pull	•	48/Case
	LBE308N	Twist-Turn	•	48/Case
900 mL	AS307N ©	Push-Pull		48/Case
	AS309N ©	Twist-Turn		48/Case
	LBE307N	Push-Pull	•	48/Case
	LBE309N	Twist-Turn	•	48/Case

### Extension Tube for Leg Bags

(HCPCS code A4331)

- Kink-resistant tubing safeguards against fluid buildup
- Flexible tubing permits leg bag to be positioned on the thigh or calf
- Sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	QTY/UNIT
AS311	18 in. (45.72 cm) Extension Tube	50/Case

### Pediatric Urine Collector

(HCPCS code A4328)

- Adhesive surface folds to form a spill-proof specimen transport container
- Hypoallergenic feature makes it safe to use with most pediatric patients
- Visible graduation marks designed for ease of use
- Single-use disposable
- Sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	QTY/UNIT
AS409 ©	200 mL Pediatric Urine Bag	50/Box

### Irrigation Trays and Kits

AMSure® Irrigation Kits and Trays are designed for simplicity and efficiency of care. The trays are available sterile or non-sterile to suit your clinical need. (HCPCS code A4320)

#### Bulb Syringe

- Flexible soft bulb promotes one-handed administration
- Designed to stand on the bulb end to help minimize contamination
- Patient label provides convenient patient documentation
- DEHP-free, BPA-free, not manufactured with natural rubber latex
- Includes syringe tip protector



AS121

#### Thumb Control Ring Syringe

- Thumb control ring top allows for maximum control and one-handed administration
- Precision crafted O-ring plunger ensures a positive seal
- Re-sealable pole bag designed for safe and convenient storage
- Patient label provides convenient patient documentation
- DEHP-free, BPA-free, not manufactured with natural rubber latex
- Includes syringe tip adapter and protector cap



AS126

Non-sterile irrigation kits include: 500 mL graduated container, patient ID label and tip protector

Sterile irrigation trays include: 1000 mL outer tray, 500 mL graduated container, waterproof underpad, alcohol prep pad, patient ID label and tip protector

#### Non-sterile Irrigation Trays

REF	THUMB CONTROL RING SYRINGE	BULB SYRINGE	RESEALABLE POLE BAG	QTY/UNIT
AS121		•		30/Case
AS126	•			30/Case
AS127	•		•	30/Case

#### Sterile Irrigation Trays

REF	THUMB CONTROL RING SYRINGE	BULB SYRINGE	QTY/UNIT
AS130 ©		•	20/Case
AS136 ©	•		20/Case



AS136



AS307N



AS121



AS126



AS136



## Ancillary & Patient Care Products

### Syringes

(HCPCS code A4322)

#### Bulb Syringe

- Flexible soft bulb promotes one-handed administration
- Designed to stand on the bulb end to help minimize contamination
- DEHP-free, BPA-free, not manufactured with natural rubber latex
- Includes syringe tip protector

#### Piston Thumb Control Ring Syringe

- Piston thumb control ring top allows for maximum control and one-handed administration
- Precision crafted O-ring plunger ensures a positive seal
- Re-sealable pole bag designed for safe and convenient storage
- DEHP-free, BPA-free, not manufactured with natural rubber latex
- Includes syringe luer tip adapter and protector cap

#### Flat Top Syringe

- Visible dosing marks designed for ease of use
- Precision crafted O-ring plunger ensures a positive seal
- Visible dosing marks designed for ease of use
- DEHP-free, BPA-free, not manufactured with natural rubber latex
- Available in sterile or non-sterile packaging
- Includes syringe luer tip adapter and protector cap

### AMSure® Specimen Containers

AMSure 4 oz., leak-resistant specimen containers are designed for the collection and safe transport of patient samples.  
(HCPCS code P9615)

- DEHP-free, BPA-free, not manufactured with natural rubber latex

REF	PACKAGING	STERILITY	LABELS	QTY/UNIT
<b>AS340</b> ©	Poly Pouch	Interior	Patient Information	100/Case
<b>AS341</b> ©	Bulk	Non-Sterile	Tamper Evident	100/Case
<b>AS343</b>	Bulk	Non-Sterile	None	500/Case
<b>AS345</b> ©	Form-Fill & Seal	Fully Sterile	None	100/Case



AS340

AS010

AS016

AS115

TYPE	REF	PACKAGING	STERILITY	TIP PROTECTOR CAP	SMALL TUBE ADAPTER	QTY/UNIT
Bulb	<b>AS011</b> ©	Poly Pouch	Sterile (EO)	•		50/Case
	<b>AS010</b>	Bulk	Non-Sterile			50/Case
Piston Thumb Control Ring	<b>AS015</b> ©	Poly Pouch	Sterile (EO)	•	•	50/Case
	<b>AS016</b>	Resealable Hanger Bag	Non-Sterile	•	•	30/Case
Piston Flat Top	<b>AS115</b> ©	Poly Pouch	Sterile (EO)	•	•	50/Case
	<b>AS116</b> ©	Resealable Pole Bag	Non-Sterile	•	•	30/Case





## Respiratory Masks and Nasal Cannulas



### AMSure® ROLLINS7™

AMSure Rollins7 is a replacement mask for virtually every type of oxygen mask on the market with its ability to provide flow ranges from 5-15 liters per minute. Rollins7 is designed with the ability to administer aerosol and nebulizer medication without interrupting oxygen therapy helping improve patient safety, patient comfort and clinical efficiency.

- Medication delivery without sacrificing FiO2. Efficient FiO2 escalation and de-escalation management
- Oxygen tubing provides versatility by connecting to flow meters with or without a Christmas tree nipple nut adapter
- Access for hydration, oral suction or oral medicine
- Non-sterile, not manufactured with natural rubber latex

REF	DESCRIPTION	QTY/UNIT
OM001 ©	Mask, Adapter, 7 ft. Star-shaped Lumen Oxygen Tubing with Label	50/Case

### Oxygen Mask

AMSure Oxygen Mask provides medium-flow oxygen delivery. The oxygen mask features smooth edges and an adjustable nose clip designed to improve patient comfort while reducing irritation.

- 7 ft. kink-resistant, star-shaped lumen tubing ensures uninterrupted flow of oxygen
- Non-sterile, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	QTY/UNIT
AS74010 ©	Adult Oxygen Mask, 7 ft. Star-shaped Lumen Oxygen Tubing	50/Case



## Enhancing Oxygen Delivery

AMSure Respiratory Therapy Systems are designed to deliver optimal FiO2 levels as well as keep upper airways unobstructed and clear to help improve quality of care. The product portfolio includes:

- Respiratory Masks and Nasal Cannulas..... 42
- Nebulizer Kits ..... 45
- Accessories..... 46
- Aerosol Drainage Bag ..... 46
- Heat and Moisture Exchangers ..... 46
- Tracheostomy Clean and Care Trays ..... 47
- Sterile Water ..... 48
- Suction Catheters, Trays and Kits ..... 49
- AMSure Suction Tubing with connectors ..... 51
- Suction Yankauers and Kits ..... 52



## Respiratory Masks and Nasal Cannulas

### High Flow Nasal Oxygen Cannula

Amsino's AMSure® high flow oxygen nasal cannula is a patient interface for delivery of humidified respiratory gases.

- Deliver flows of up to 60 lpm
- Rotatable silicone clip for optimizing patient comfort and fit
- Flexible EVA tube with detachable and slide-in pliable nasal prongs for connection to device on either side of patient
- Universal connector to be applied with most high flow nasal cannula devices
- Single patient use, not manufactured with natural rubber latex

REF	DESCRIPTION	FEATURES	QTY/UNIT
<b>AS75280</b> ©	High Flow Nasal Oxygen Cannula	Small, Adult	10/Case
<b>AS75285</b>	High Flow Nasal Oxygen Cannula	Medium, Adult	10/Case
<b>AS75282</b> ©	High Flow Nasal Oxygen Cannula	Large, Adult	10/Case



AS75280

### Nasal Oxygen Cannula

AMSure Nasal Oxygen Cannulas provide low-to-medium flow oxygen delivery and are designed with soft and pliable prongs which are ideal for patients who require oxygen therapy.

- Over the ear design allows for complete freedom of movement
- 7 ft. kink-resistant, star-shaped lumen tubing ensures uninterrupted flow of oxygen
- Non-sterile, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	FEATURES	QTY/UNIT
<b>AS75080</b> ©	Standard Nasal Oxygen Cannula	Flared Prongs, Adult	50/Case
<b>AS75085</b> ©	Standard Nasal Oxygen Cannula	Non-Flared Prongs, Adult	50/Case



AS75080

### Cushion Nasal Oxygen Cannula

AMSure Cushion Nasal Oxygen Cannulas with soft foam ear covers provide low-to-medium flow oxygen delivery. Soft, pliable prongs are ideal for patients who require long-term oxygen therapy.

- Over the ear design allows for complete freedom of movement
- 7 ft. kink-resistant, star-shaped lumen tubing ensures uninterrupted flow of oxygen
- Non-sterile, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	FEATURES	QTY/UNIT
<b>AS75181</b> ©	Soft Nasal Oxygen Cannula	Curved Prongs with Cushion, Adult	25/Case




### Soft Nasal Oxygen Cannula

AMSure® Soft Nasal Oxygen Cannulas provide low-to-medium flow oxygen delivery and are designed with soft and pliable prongs which are ideal for patients who require long-term oxygen therapy.

- Over the ear design allows for complete freedom of movement
- 7 ft. kink-resistant, star-shaped lumen tubing ensures uninterrupted flow of oxygen
- Non-sterile, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	FEATURES	QTY/UNIT
<b>AS75180</b> ©	Soft Nasal Oxygen Cannula	Curved Prongs, Adult	50/Case
<b>AS75190</b> ©	Soft Nasal Oxygen Cannula	Curved Prongs, Pediatric	50/Case
<b>AS75199</b> ©	Soft Nasal Oxygen Cannula	Curved Prongs, Infant	50/Case



AS75180

### Non-Rebreather Oxygen Mask

AMSure Non-rebreather Oxygen Masks provide high-flow oxygen delivery. The non-rebreather mask features smooth edges and an adjustable nose clip designed to improve patient comfort while reducing irritation.

- Soft and comfortable elastic strap reduces facial pressure while providing a tight seal
- 7 ft. kink-resistant, star-shaped lumen tubing ensures uninterrupted flow of oxygen
- Non-sterile, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	QTY/UNIT
<b>AS75010</b> ©	Adult Non-rebreather Oxygen Mask, 7 ft. Star-shaped Lumen Oxygen Tubing	50/Case
<b>AS75020</b> ©	Pediatric Non-rebreather Oxygen Mask, 7 ft. Star-shaped Lumen Oxygen Tubing	50/Case



AS75010

### Aerosol Mask

AMSure Aerosol Mask for pediatric patients fits any standard 22 mm tubing or nebulizer. The aerosol mask features smooth edges and an adjustable nose clip designed to improve patient comfort while reducing irritation.

- Soft and comfortable elastic strap reduces facial pressure
- Non-sterile, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	QTY/UNIT
<b>AS75040</b> ©	Adult Aerosol Mask	50/Case
<b>AS75050</b> ©	Pediatric Aerosol Mask	50/Case



AS75050

## Respiratory Masks and Nasal Cannulas

### Tracheostomy Mask

AMSure® Tracheostomy Mask features a soft cuff to improve patient comfort. The 22 mm connector fits standard corrugated tubing for humidified oxygen and nebulizer treatments.

- Soft cuff reduces pressure and irritation to patient skin
- Adapter swivels and rotates 360 degrees, ensuring optimum position for attachments
- Non-sterile, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	QTY/UNIT
AS75060 ©	Adult Tracheostomy Mask	50/Case



### Nebulizer Kits

AMSure Nebulizer Kits are designed for ease of use and optimal particle size delivery. The nebulizer kits feature a 6 mL medication cup, 7 ft. tubing and either a T-mouthpiece with corrugated tubing or curved mouthpiece.


- Anti-spill medication container lid helps to prevent leakage
- Visible dosing marks designed for ease of use
- 7 ft. kink-resistant, star-shaped lumen tubing ensures uninterrupted flow of oxygen
- Non-sterile, BPA-free, not manufactured with natural rubber latex

#### T-Mouthpiece

- Disposable nebulizers deliver effective treatment to a wide range of patients

#### Curved Mouthpiece

- Disposable and reusable nebulizers delivers effective treatment to a wide range of patients



AS78010



AS78015

REF	DESCRIPTION	FEATURES	QTY/UNIT
AS78010 ©	Nebulizer, 5 Piece, 7 ft. Star-shaped Lumen Oxygen Tubing	T-mouthpiece	50/Case
AS78015 ©	Nebulizer, 3 Piece, 7 ft. Star-shaped Lumen Oxygen Tubing	Curved mouthpiece	50/Case
AS78020 ©	Nebulizer, 3 Piece, 7 ft. Star-shaped Lumen Oxygen Tubing	Aerosol mask	50/Case
AS78025 ©	Nebulizer, 3 Piece, 7 ft. Star-shaped Lumen Oxygen Tubing	Curved mouthpiece, Reusable	50/Case

## Accessories

### Oxygen Tubing

7 ft. Star-shaped Lumen Tubing with Standard 22 mm Connector.

- 7 ft. kink-resistant, star-shaped lumen tubing ensures uninterrupted flow of oxygen
- Non-sterile, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	FEATURES	QTY/UNIT
AS76007 ©	7 ft. Oxygen Tubing	Star-shaped Lumen	50/Case



### Respiratory Adapters

AMSure® Respiratory Adapters are available in 15 mm I.D. and 22 mm I.D.

- Non-sterile, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	FEATURES	QTY/UNIT
AS730100	Universal Cuff Adapter	22 mm ID Connector	50/Case
AS730200	Multi Adapter	15 mm ID / 22 mm OD Connector	50/Case



AS730100

AS730200

## Aerosol Drainage Bag

AMSure Aerosol Drainage System collects moisture accumulating from continuous in-line aerosol therapy. The system includes a 750 mL drainage bag, Y adapter, self-locking plastic chain and accepts 22 mm I.D. corrugated tubing.

- Pre-connected self-locking plastic chain secures the system to a stable object
- Convenient drain port allows for easy one-handed drainage
- Non-sterile, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	QTY/UNIT
AS75035	750 mL Aerosol Drainage Bag	50/Case





## Heat and Moisture Exchangers



HME001

Amsino is offering Passive Heat and Moisture Exchanger (HME) for Tracheostomy Patients. HME is required when a patient can no longer inspire and expire through the nose and mouth. The nose and mouth provide natural warmth, moisture and filtration for the air we breathe. Having a tracheostomy tube, however, by-passes these mechanisms so humidification must be provided to keep secretions thin and to avoid mucus plugs. In order to keep the environment at an optimal humidity level, Heat Moisture Exchangers are used.

REF	DESCRIPTION	QTY/UNIT
<b>HME001</b> ©	Heat and Moisture Exchangers, Non-Sterile	50/Case
<b>HME002</b> ©	Heat and Moisture Exchangers, Sterile	50/Case

## AMSure® Sterile Water

All products are DEHP-Free, PVC-Free, not manufactured with natural rubber latex, preservative-free, and Non-Pyrogenic.



AS1065

AS1000

### Sterile Water for Inhalation (Bottles)

- 500 mL and 750 mL, and 1000 mL bottles Sterile Water, USP
- Used in nebulizers and humidifiers

REF	REF (Outside US only)	DESCRIPTION	QTY/UNIT
<b>AS1565</b>	<b>AS1565E</b> ©	500 mL Sterile Water for Inhalation, USP	12/Case
<b>AS1765</b>	<b>AS1765E</b> ©	750 mL Sterile Water for Inhalation, USP	12/Case
<b>AS1065</b>	<b>AS1065E</b> ©	1000 mL Sterile Water for Inhalation, USP	12/Case
<b>AS1066</b>	—	500 mL Sterile Water for Inhalation with non-sterile cap	12/Case
<b>AS1566</b>	—	1000 mL Sterile Water for Inhalation with non-sterile cap	12/Case
<b>AS1000</b>	—	1000 mL Sterile Water for Inhalation, Hanging Spikeable Bottle	12/Case

### Sterile Water for Inhalation (Vials)

- Available in two sizes: 3 mL and 5 mL, USP
- Easy open twist-off cap is simple to use and allows dispensing by the drop or milliliter
- Used for suctioning and respiratory inhalation



ASR0035

REF	REF (Outside US only)	DESCRIPTION	QTY/UNIT
<b>ASR0035</b>	<b>ASR0035E</b> ©	3 mL Sterile Water Inhalation, USP	1000/Case
<b>ASR0055</b>	<b>ASR0055E</b> ©	5 mL Sterile Water Inhalation, USP	1000/Case

### Sterile Water Prefilled Humidifiers

- 350 mL and 550 mL bottles Sterile Water, USP for inhalation therapy and humidification
- For respiratory therapy, adapters included with each bottle
- Reservoir bottles have built-in oxygen tubing connectors
- Micro-diffuser produces smaller bubbles and greater surface agitation for quiet operation
- Audible alarm alerts clinician of flow restriction or occlusion of oxygen tubing



AS0352

AS0552

REF	REF (Outside US only)	DESCRIPTION	QTY/UNIT
<b>AS0352</b>	<b>AS0352E</b> ©	350 mL Sterile Water for Inhalation Prefilled Humidifier, USP	20/Case
<b>AS0552</b>	<b>AS0552E</b> ©	550 mL Sterile Water for Inhalation Prefilled Humidifier, USP	12/Case

## Tracheostomy Clean and Care Trays



AS861

AMSure® Tracheostomy Clean and Care Trays feature the necessary components to perform cleaning in and around the stoma. The trays provide clinicians the versatility to perform care based on the patient's individual needs.

- Sterile, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	FEATURES	QTY/UNIT
<b>AS860</b>	Cleaning Brush, Vinyl Powder-free Gloves, (2) Cotton Tip Applicators with Plastic Shaft, Twill Tape, Tracheostomy Dressing and Waterproof Drape	Tray, (2) 4" x 4" Gauze, (2) Pipe Cleaners, Pop-up Solution Cup	20/Case
<b>AS861</b>	Cleaning Brush, Vinyl Powder-free Gloves, (2) Cotton Tip Applicators with Plastic Shaft, Twill Tape, Tracheostomy Dressing and Waterproof Drape	2 Compartment Tray, Removable Inner Tray, (4) 4" x 4" Gauze, (4) Pipe Cleaners	20/Case

# Suction Catheters, Trays and Kits

## Suction Catheters

AMSure® Suction Catheters improve patient care by safely removing respiratory secretions from the upper airway. The suction catheters are available in both straight and coiled packs in various French sizes.

- **Opposing staggered eyelets reduce the potential for blocked suction**
- **Whistle cap port enhances controlled suction**
- **Sterile, BPA-free, not manufactured with natural rubber latex**



AS361

## Suction Catheter Kits

AMSure Suction Catheters improve patient care by safely removing respiratory secretions from the upper airway. The suction catheters are available in both straight and coiled packs in various French sizes.

- **Uniquely designed two-eyed whistle tip offers safe and effective suctioning**
- **Sterile, BPA-free, not manufactured with natural rubber latex**



AS371

REF	DESCRIPTION	FEATURES	QTY/UNIT
AS361 ©	6 Fr, 17 ½ in.	Straight Pack	50/Case
AS361C ©	6 Fr, 17 ½ in.	Coiled Pack	50/Case
AS362 ©	8 Fr, 17 ½ in.	Straight Pack	50/Case
AS362C ©	8 Fr, 17 ½ in.	Coiled Pack	50/Case
AS363 ©	10 Fr, 23 in.	Straight Pack	50/Case
AS363C ©	10 Fr, 23 in.	Coiled Pack	50/Case
AS364 ©	12 Fr, 23 in.	Straight Pack	50/Case
AS364C ©	12 Fr, 23 in.	Coiled Pack	50/Case
AS365 ©	14 Fr, 23 in.	Straight Pack	50/Case
AS365C ©	14 Fr, 23 in.	Coiled Pack	50/Case
AS366 ©	16 Fr, 23 in.	Straight Pack	50/Case
AS366C ©	16 Fr, 23 in.	Coiled Pack	50/Case
AS367 ©	18 Fr, 23 in.	Straight Pack	50/Case
AS367C ©	18 Fr, 23 in.	Coiled Pack	50/Case

REF	DESCRIPTION	FEATURES	QTY/UNIT
AS371 ©	6 Fr, 17 ½ in.	Whistle Tip	50/Case
AS381 ©	6 Fr, 17 ½ in.	Whistle Tip	50/Case
AS372 ©	8 Fr, 17 ½ in.	Whistle Tip	50/Case
AS382 ©	8 Fr, 17 ½ in.	Whistle Tips	50/Case
AS373 ©	10 Fr, 23 in.	Whistle Tip	50/Case
AS383 ©	10 Fr, 23 in.	Whistle Tip	50/Case
AS374 ©	12 Fr, 23 in.	Whistle Tip	50/Case
AS384 ©	12 Fr, 23 in.	Whistle Tip	50/Case
AS375 ©	14 Fr, 23 in.	Whistle Tip	50/Case
AS385 ©	14 Fr, 23 in.	Whistle Tip	50/Case
AS376 ©	16 Fr, 23 in.	Whistle Tip	50/Case
AS386 ©	16 Fr, 23 in.	Whistle Tip	50/Case
AS387 ©	18 Fr, 23 in.	Whistle Tip	50/Case

## Mini Suction Catheter Trays

The AMSure® Mini Suction Catheter Trays' space-saving design offers everything needed for effective procedures. The mini tray serves as a basin to hold sterile water or saline.



SCT10


- **Opposing staggered eyelets reduce the potential for blocked suction**
- **DeLee Tip's smooth and rounded tip minimizes trauma for pediatric patients**
- **Whistle tip enhances controlled suction**
- **Sterile, BPA-free, not manufactured with natural rubber latex**

REF	DESCRIPTION	FEATURES	QTY/UNIT
SCT08	8 Fr, 17 ½ in. Suction Catheter, Tray and Vinyl Gloves	DeLee Tip	100/Case
SCT10	10 Fr, 23 in. Suction Catheter, Tray and Vinyl Gloves	DeLee Tip	100/Case
SCT12	12 Fr, 23 in. Suction Catheter, Tray and Vinyl Gloves	Whistle Tip	100/Case
SCT14	14 Fr, 23 in. Suction Catheter, Tray and Vinyl Gloves	Whistle Tip	100/Case
SCT16	16 Fr, 23 in. Suction Catheter, Tray and Vinyl Gloves	Whistle Tip	100/Case

## Closed Suction Catheter Trays

24-hour Adult Catheter from 10 FR to 14 FR Endotracheal Catheter Kits and Tracheostomy Catheter Kits with MDI Port Adapter.

- **Closed system to prevent cross infections in facility**
- **Clean easily with needle free port**
- **Avoid interruptions of oxygen supply due to suction**



REF	DESCRIPTION	FLEX CONNECTOR	QTY/UNIT
Endotracheal Kits with MDI Port Adapter			
VS363RC124 ©	10 Fr, 21.3 in. Endotracheal Catheter, Double Swivel Elbow	15 mm x 22 mm	20/Case
VS364RC124 ©	12 Fr, 21.3 in. Endotracheal Catheter, Double Swivel Elbow	15 mm x 22 mm	20/Case
VS365RC124 ©	14 Fr, 21.3 in. Endotracheal Catheter, Double Swivel Elbow	15 mm x 22 mm	20/Case
VS363TP124 ©	10 Fr, 21.3 in. Endotracheal Catheter, T-Piece	22 mm x 22 mm	20/Case
VS364TP124 ©	12 Fr, 21.3 in. Endotracheal Catheter, T-Piece	22 mm x 22 mm	20/Case
VS365TP124 ©	14 Fr, 21.3 in. Endotracheal Catheter, T-Piece	22 mm x 22 mm	20/Case
Tracheostomy Kits with MDI Port Adapter			
VS363TC124 ©	10 Fr, 12 in. Tracheostomy Catheter, Double Swivel Elbow	15 mm x 22 mm	20/Case
VS364TC124 ©	12 Fr, 12 in. Tracheostomy Catheter, Double Swivel Elbow	15 mm x 22 mm	20/Case
VS365TC124 ©	14 Fr, 12 in. Tracheostomy Catheter, Double Swivel Elbow	15 mm x 22 mm	20/Case
VS363TTP124 ©	10 Fr, 12 in. Tracheostomy Catheter, T-Piece	22 mm x 22 mm	20/Case
VS364TTP124 ©	12 Fr, 12 in. Tracheostomy Catheter, T-Piece	22 mm x 22 mm	20/Case
VS365TTP124 ©	14 Fr, 12 in. Tracheostomy Catheter, T-Piece	22 mm x 22 mm	20/Case



## Suction Tubing with Connectors

AMSure® Suction Tubing, with molded female connectors and male adapters, provide secure connections to support uninterrupted procedures. Flexible, low-coil memory allows for easy draping.

- Sterile and non-sterile options, BPA-free, not manufactured with natural rubber latex



AS825

DIAMETER (ID)	LENGTH	STERILE	NON-STERILE
3/16 in. (5 mm) Green Connectors	18 in. (46 cm)	<b>AS820</b> © (50/Case)	<b>AS820B</b> © (100/Case)
	6 ft. (1.8 m)	<b>AS821</b> © (50/Case)	<b>AS821B</b> © (50/Case)
	10 ft. (3.0 m)	<b>AS822</b> © (50/Case)	<b>AS822B</b> © (50/Case)
	12 ft. (3.7 m)	<b>AS823</b> © (20/Case)	<b>AS823B</b> © (50/Case)
	20 ft. (6.1 m)	<b>AS824</b> © (20/Case)	<b>AS824B</b> © (50/Case)
1/4 in. (6 mm) Light Blue Connectors	6 ft. (1.8 m)	<b>AS825</b> © (50/Case)	<b>AS825B</b> © (50/Case)
	10 ft. (3.0 m)	<b>AS826</b> © (50/Case)	<b>AS826B</b> © (50/Case)
	10 ft. (6 mm × 3.0 m) with Pinch Clamps	<b>AS826C</b> © (50/Case)	—
	12 ft. (3.7 m)	<b>AS827</b> © (20/Case)	<b>AS827B</b> © (50/Case)
	20 ft. (6.1 m)	<b>AS828</b> © (20/Case)	<b>AS828B</b> © (50/Case)
	100 ft. (30.0 m) with Connectors every 6 ft. (1.8 m)	—	<b>AS829B-01</b> © (1Roll/Case)
9/32 in. (7 mm) Dark Blue Connectors	6 ft. (1.8 m)	<b>AS825M</b> © (50/Case)	<b>AS825BM</b> © (50/Case)
	10 ft. (3.0 m)	<b>AS826M</b> © (50/Case)	<b>AS826BM</b> © (50/Case)
	12 ft. (3.7 m)	<b>AS827M</b> © (20/Case)	<b>AS827BM</b> © (50/Case)
	20 ft. (6.1 m)	<b>AS828M</b> © (20/Case)	<b>AS828BM</b> © (50/Case)

## Suction Yankauers and Kits

AMSure® Suction Yankauers are designed with a slip-resistant, lightweight handle that helps provide optimal control and maneuverability.

- Translucent pathway allows for fluid visualization
- Smooth tips are designed with side eyes to minimize the chance of tissue grab and damage
- Sterile, disposable, DEHP-free, BPA-free, not manufactured with natural rubber latex



AS835

AS832

### Yankauers

DESCRIPTION			STERILE	NON-STERILE
Rigid	Regular Beveled Tip	Non-Vented	<b>AS830</b> © 50/Case	<b>AS830B</b> 500/Case
Rigid	Regular Beveled Tip	Vented	<b>AS831</b> © 50/Case	<b>AS831B</b> 500/Case
Rigid	Bulbous Tip	Non-Vented	<b>AS832</b> © 50/Case	<b>AS832B</b> 500/Case
Rigid	Bulbous Tip	Vented	<b>AS833</b> © 50/Case	<b>AS833B</b> 500/Case
Flexible	Regular Beveled Tip	Non-Vented	<b>AS834</b> © 50/Case	<b>AS834B</b> 500/Case
Flexible	Regular Beveled Tip	Vented	<b>AS835</b> © 50/Case	<b>AS835B</b> 500/Case
Flexible	Bulbous Tip	Non-Vented	<b>AS836</b> © 50/Case	<b>AS836B</b> 500/Case
Flexible	Bulbous Tip	Vented	<b>AS837</b> © 50/Case	<b>AS837B</b> 500/Case

### Yankauer & Tubing Kits

REF	DESCRIPTION	QTY/UNIT
<b>AS825P</b> ©	Rigid Yankauer, Bulbous Tip, Non-Vented, pre-attached to 1/4 in. (6 mm) × 6 ft. (1.8 m) Suction Tubing, Sterile	50/Case
<b>AS827P*</b> ©	Rigid Yankauer, Non-Vented with Bulbous Tip, pre-attached to 1/4 in. (6 mm) × 12 ft. (3.7 m) Suction Tubing, Sterile	50/Case

\*Made to order

# AMSure®

## STERILE WATER AND SALINE PRODUCTS



## Inhalation, Irrigation, and Injection Solutions

AMSure Respiratory Therapy Systems are designed to deliver optimal FiO2 levels as well as keep upper airways unobstructed and clear to help improve quality of care. The product portfolio includes:

- Sterile Water for Inhalation, USP ..... 54
- Sterile Water Prefilled Humidifiers, USP ..... 54
- Nebulizer Kits ..... 45
- Saline Products ..... 55
- Device and Wound Irrigation ..... 56

## AMSure® Sterile Water

All products are DEHP-Free, PVC-Free, not manufactured with natural rubber latex, preservative-free, and Non-Pyrogenic.



AS1065

AS1000

### Sterile Water for Inhalation (Bottles)

- 500 mL and 750 mL, and 1000 mL bottles Sterile Water, USP
- Used in nebulizers and humidifiers

REF	REF (Outside US only)	DESCRIPTION	QTY/UNIT
AS1565	AS1565E ©	500 mL Sterile Water for Inhalation, USP	12/Case
AS1765	AS1765E ©	750 mL Sterile Water for Inhalation, USP	12/Case
AS1065	AS1065E ©	1000 mL Sterile Water for Inhalation, USP	12/Case
AS1066	—	500 mL Sterile Water for Inhalation with non-sterile cap	12/Case
AS1566	—	1000 mL Sterile Water for Inhalation with non-sterile cap	12/Case
AS1000	—	1000 mL Sterile Water for Inhalation, Hanging Spikeable Bottle	12/Case

### Sterile Water for Inhalation (Vials)

- Available in two sizes: 3 mL and 5 mL, USP
- Easy open twist-off cap is simple to use and allows dispensing by the drop or milliliter
- Used for suctioning and respiratory inhalation



ASR0035

REF	REF (Outside US only)	DESCRIPTION	QTY/UNIT
ASR0035	ASR0035E ©	3 mL Sterile Water Inhalation, USP	1000/Case
ASR0055	ASR0055E ©	5 mL Sterile Water Inhalation, USP	1000/Case

### Sterile Water Prefilled Humidifiers

- 350 mL and 550 mL bottles Sterile Water, USP for inhalation therapy and humidification
- For respiratory therapy, adapters included with each bottle
- Reservoir bottles have built-in oxygen tubing connectors
- Micro-diffuser produces smaller bubbles and greater surface agitation for quiet operation
- Audible alarm alerts clinician of flow restriction or occlusion of oxygen tubing



AS0352

AS0552

REF	REF (Outside US only)	DESCRIPTION	QTY/UNIT
AS0352	AS0352E ©	350 mL Sterile Water for Inhalation Prefilled Humidifier, USP	20/Case
AS0552	AS0552E ©	550 mL Sterile Water for Inhalation Prefilled Humidifier, USP	12/Case



# Saline Products

## AMSafe® Prefilled Saline Flush Syringe

- 3 mL, 5 mL, 10 mL and 20 mL syringes, 0.9% Sodium Chloride Injection, USP
- Sterile fluid path (dust cover) and sterile field configurations
- Complies with PICC and INS standards
- Compatible with standard ISO Luer connectors
- FDA cleared
- Sterile, non-pyrogenic, DEHP-free, PVC-free, not manufactured with natural rubber latex



## 12 mL and 20 mL Syringes

### Standard Dust Cover Flush Syringes (Sterile Fluid Pathway)

REF	DESCRIPTION	QTY/UNIT
IVF1203	3 mL USP Sodium Chloride Injection Fill in 12 mL Syringe	180/Case (6 Boxes – 30 per Box)
IVF1205	5 mL USP Sodium Chloride Injection Fill in 12 mL Syringe	180/Case (6 Boxes – 30 per Box)
IVF1210	10 mL USP Sodium Chloride Injection Fill in 12 mL Syringe	180/Case (6 Boxes – 30 per Box)
IVF2020	20 mL USP Sodium Chloride Injection Fill in 20 mL Syringe	120/Case (6 Boxes – 20 per Box)

### Sterile Field Ready Flush Syringes

REF	DESCRIPTION	QTY/UNIT
IVF1203SS	3 mL USP Sodium Chloride Injection Fill in 12 mL Syringe	180/Case (6 Boxes – 30 per Box)
IVF1205SS	5 mL USP Sodium Chloride Injection Fill in 12 mL Syringe	180/Case (6 Boxes – 30 per Box)
IVF1210SS	10 mL USP Sodium Chloride Injection Fill in 12 mL Syringe	180/Case (6 Boxes – 30 per Box)
IVF2020SS	20 mL USP Sodium Chloride Injection Fill in 20 mL Syringe	120/Case (6 Boxes – 20 per Box)



## 10 mL Syringes

### Standard Dust Cover Flush Syringes (Sterile Fluid Pathway)

REF	DESCRIPTION	QTY/UNIT
IVF1003	3 mL USP Sodium chloride Injection Fill in 10 mL Syringe	180/Case (6 Boxes – 30 per Box)
IVF1005	5 mL USP Sodium chloride Injection Fill in 10 mL Syringe	180/Case (6 Boxes – 30 per Box)
IVF1010	10 mL USP Sodium chloride Injection Fill in 10 mL Syringe	240/Case (8 Boxes – 30 per Box)
IVF1010T	10 mL USP Sodium chloride Injection Fill in 10 mL Syringe	180/Case (6 Boxes – 30 per Box)
IVF1010TM	10 mL USP Sodium chloride Injection Fill in 10 mL Syringe	240/Case (8 Boxes – 30 per Box)

## AMSure® Sodium Chloride Inhalation Solution (vials)

- 3 mL, 5 mL and 15 mL Sodium Chloride Inhalation Solution (0.9%), USP
- Single-use vials with graduations for accurate dispensing
- Easy open twist-off cap is simple to use and allows dispensing by the drop or milliliter
- Used for suctioning and respiratory inhalation

REF	REF (Outside US only)	DESCRIPTION	QTY/UNIT
ASR0039	ASR0039E	3 mL Sodium Chloride Inhalation Solution (0.9%) USP	1000/Case
ASR0059	ASR0059E	5 mL Sodium Chloride Inhalation Solution (0.9%) USP	1000/Case
ASR0159	ASR0159E	15 mL Sodium Chloride Inhalation Solution (0.9%) USP	144/Case



ASR0039

ASR0059

ASR0159

## AMSure Device and Wound Irrigation

- 100 mL sterile water or saline bottle with screw cap and inner seal
- Easy open twist screw top container
- Used for device and wound irrigation



AS399

AS398

REF	REF (Outside US only)	DESCRIPTION	QTY/UNIT
AS399	AS399E	100 mL Sterile Water for Device and Wound Irrigation	48/Case
AS398	AS398E	100 mL Sterile Saline 0.9% for Device and Wound Irrigation	48/Case



### Suction Tubing with Connectors

AMSure® Suction Tubing, with molded female connectors and male adapters, provide secure connections to support uninterrupted procedures. Flexible, low-coil memory allows for easy draping.

- Sterile and non-sterile options, BPA-free, not manufactured with natural rubber latex



AS825

DIAMETER (ID)	LENGTH	STERILE	NON-STERILE
3/16 in. (5 mm) Green Connectors	18 in. (46 cm)	<b>AS820</b> © (50/Case)	<b>AS820B</b> © (100/Case)
	6 ft. (1.8 m)	<b>AS821</b> © (50/Case)	<b>AS821B</b> © (50/Case)
	10 ft. (3.0 m)	<b>AS822</b> © (50/Case)	<b>AS822B</b> © (50/Case)
	12 ft. (3.7 m)	<b>AS823</b> © (20/Case)	<b>AS823B</b> © (50/Case)
	20 ft. (6.1 m)	<b>AS824</b> © (20/Case)	<b>AS824B</b> © (50/Case)
1/4 in. (6 mm) Light Blue Connectors	6 ft. (1.8 m)	<b>AS825</b> © (50/Case)	<b>AS825B</b> © (50/Case)
	10 ft. (3.0 m)	<b>AS826</b> © (50/Case)	<b>AS826B</b> © (50/Case)
	10 ft. (6 mm × 3.0 m) with Pinch Clamps	<b>AS826C</b> © (50/Case)	—
	12 ft. (3.7 m)	<b>AS827</b> © (20/Case)	<b>AS827B</b> © (50/Case)
	20 ft. (6.1 m)	<b>AS828</b> © (20/Case)	<b>AS828B</b> © (50/Case)
	100 ft. (30.0 m) with Connectors every 6 ft. (1.8 m)	—	<b>AS829B-01</b> © (1Roll/Case)
9/32 in. (7 mm) Dark Blue Connectors	6 ft. (1.8 m)	<b>AS825M</b> © (50/Case)	<b>AS825BM</b> © (50/Case)
	10 ft. (3.0 m)	<b>AS826M</b> © (50/Case)	<b>AS826BM</b> © (50/Case)
	12 ft. (3.7 m)	<b>AS827M</b> © (20/Case)	<b>AS827BM</b> © (50/Case)
	20 ft. (6.1 m)	<b>AS828M</b> © (20/Case)	<b>AS828BM</b> © (50/Case)

### Everyday Patient Care Solutions

AMSure Surgical and Nursing Systems offer a variety of frequently used products for patient care. These essentials include:

- Suction Tubing ..... 58
- Suction Yankauers and Kits ..... 59
- Suction Device Holster System ..... 59
- Specimen Containers ..... 60
- Vaginal Speculums ..... 60
- Disposable Umbilical Cord Clamp ..... 60
- Ear/Ulcer Syringes ..... 61
- Cleansing Enema Sets ..... 61
- Connectors ..... 62



## Suction Yankauers and Kits

AMSure® Suction Yankauers are designed with a slip-resistant, lightweight handle that helps provide optimal control and maneuverability.

- Translucent pathway allows for fluid visualization
- Smooth tips are designed with side eyes to minimize the chance of tissue grab and damage
- Sterile, disposable, DEHP-free, BPA-free, not manufactured with natural rubber latex



AS835 AS832

### Yankauers

DESCRIPTION			STERILE	NON-STERILE
Rigid	Regular Beveled Tip	Non-Vented	<b>AS830</b> © 50/Case	<b>AS830B</b> 500/Case
Rigid	Regular Beveled Tip	Vented	<b>AS831</b> © 50/Case	<b>AS831B</b> 500/Case
Rigid	Bulbous Tip	Non-Vented	<b>AS832</b> © 50/Case	<b>AS832B</b> 500/Case
Rigid	Bulbous Tip	Vented	<b>AS833</b> © 50/Case	<b>AS833B</b> 500/Case
Flexible	Regular Beveled Tip	Non-Vented	<b>AS834</b> © 50/Case	<b>AS834B</b> 500/Case
Flexible	Regular Beveled Tip	Vented	<b>AS835</b> © 50/Case	<b>AS835B</b> 500/Case
Flexible	Bulbous Tip	Non-Vented	<b>AS836</b> © 50/Case	<b>AS836B</b> 500/Case
Flexible	Bulbous Tip	Vented	<b>AS837</b> © 50/Case	<b>AS837B</b> 500/Case

### Yankauer & Tubing Kits

REF	DESCRIPTION	QTY/UNIT
<b>AS825P</b> ©	Rigid Yankauer, Bulbous Tip, Non-Vented, pre-attached to 1/4 in. (6 mm) x 6 ft. (1.8 m) Suction Tubing, Sterile	50/Case
<b>AS827P*</b> ©	Rigid Yankauer, Non-Vented with Bulbous Tip, pre-attached to 1/4 in. (6 mm) x 12 ft. (3.7 m) Suction Tubing, Sterile	50/Case

\*Made to order

## Specimen Containers

AMSure® 4 oz., leak-resistant specimen containers are designed for the collection and safe transport of patient samples.

- DEHP-free, BPA-free, not manufactured with natural rubber latex

REF	FORM FILL SEAL	POLY POUCH	BULK	STERILE INTERIOR	TAMPER EVIDENT LABEL	PATIENT LABELS	QTY/UNIT
<b>AS340</b> ©		•		•		•	100/Case
<b>AS341</b> ©			•	•	•		100/Case
<b>AS343</b>			•				500/Case
<b>AS345</b> ©	•			•			100/Case



AS340

## Vaginal Speculums

AMSure Vaginal Speculum helps to increase patient comfort during insertion and may provide greater visualization. The speculums are available in small, medium and large sizes.

- Round-edge design helps to increase patient comfort
- Single-use disposable
- Non-sterile, not manufactured with natural rubber latex



AS032M

REF	DESCRIPTION	QTY/UNIT
<b>AS032S</b> ©	Small Vaginal Speculum	100/Case
<b>AS032M</b> ©	Medium Vaginal Speculum	100/Case
<b>AS032L</b> ©	Large Vaginal Speculum	100/Case

## Suction Device Holster System

SteriCadi™ is a sterile disposable holder to store suction devices. The intuitive design is easy to set up and use to help support better infection control outcomes.

- Suction devices can be isolated to help reduce cross-contamination and improve safety measures
- Sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	QTY/UNIT
<b>SIH0001</b> ©	SteriCadi Suction Device Holder with Clip	30/Case



## Disposable Umbilical Cord Clamp

AMSure Umbilical Cord Clamp with safety lock helps provide a secure connection.

- Safety lock is designed to decrease risk of accidental disconnection
- Grooved teeth help provide an effective grip
- Sterile, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	QTY/UNIT
<b>UCC100</b>	Umbilical Cord Clamp, sterile	50/Case





## Ear/Ulcer Syringes

AMSure® Ear/Ulcer syringes are ideal for irrigating smaller wounds and alleviating congestion. The syringes are available in several sizes to meet clinical requirements.

- Gentle tip design reduces irritation
- Ribbed bulb provides steady grip control
- Designed to stand on the bulb end to help minimize contamination
- Sterile and non-sterile options, not manufactured with natural rubber latex

REF	DESCRIPTION	FEATURES	QTY/UNIT
AS00501	1 oz. syringe	Non-sterile, Bulk	300/Case
AS00502	2 oz. syringe	Non-sterile, Bulk	200/Case
AS00503	3 oz. syringe	Non-sterile, Bulk	150/Case
AS00501S ©	1 oz. syringe	Sterile, Single Pack	50/Case
AS00502S ©	2 oz. syringe	Sterile, Single Pack	50/Case
AS00503S ©	3 oz. syringe	Sterile, Single Pack	50/Case

## Connectors

### 5-in-1 Straight Connector

AMSure® 5-in-1 Connector is a durable connector used to connect same size tubing or use as a reducer and fits five (3/16 in., 1/4 in., 5/16 in., 3/8 in., 7/16 in.) sizes of tubing.

- DEHP-free, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	FEATURES	QTY/UNIT
AS807 ©	5-in-1 Connector	Sterile	100/Case
AS807B	5-in-1 Connector	Non-sterile, Bulk	500/Case



AS807

### 6-in-1 Universal Y Connector

Lightweight, durable polypropylene plastic connectors adapt for almost any use allowing for multiple pieces of tubing to be joined together. Tapered legs fit six sizes of tubing, from 6 mm to 15 mm (1/4 in. to 9/16 in.).

- DEHP-free, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	QTY/UNIT
UYC001S	6-in-1 Universal Y Connector, Sterile	100/Case
UYC001	6-in-1 Universal Y Connector, Non-sterile	500/Case



AS330

## Cleansing Enema Sets

AMSure Cleansing Enema Bag and Bucket Sets provide the necessary components to perform a comfortable and safe procedure. The set includes 60 in. (152.4 cm) tubing with pre-lubricated tip, soap packet, waterproof drape and either a 1500 mL bag or 1500 mL bucket.

### Bag Set

- One-way bag valve prevents solution spills
- Pre-lubricated tip allows for easy and comfortable insertion
- Non-sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex

### Bucket Set

- Graduated markings designed for precision and ease of use
- Built-in handles provide convenient transport
- Pre-lubricated tip allows for easy and comfortable insertion
- Non-sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	FEATURES	QTY/UNIT
AS330	60 in. (152.4 cm) Tubing with Pre-Lubricated Tip, Slide Clamp, Soap Packet and Waterproof Drape	1500 mL Bag	50/Case
AS333	60 in. (152.4 cm) Tubing with Pre-Lubricated Tip, Slide Clamp, Soap Packet and Waterproof Drape	1500 mL Bucket	50/Case



AS333





# RECEPTAL®

## SUCTION & WASTE FLUID MANAGEMENT SYSTEMS



## Waste Fluid Management Solutions

Amsino is a leading global provider of suction fluid management systems. The breadth of Amsino's RECEPTAL suction fluid management systems enhance patient quality of care through clinical efficacy, ease-of-use, product selection and key safety features. Our product selection includes:

- RECEPTAL Outer Shells ..... 64
- RECEPTAL Liners ..... 65
- RECEPTAL Twist Liners ..... 66
- iRECEPTAL® Surgical Waste Fluid Management System ..... 67
- EZE-VAC® Suction Canister System ..... 69
- EZE-VAC 2.0 Suction Canister System ..... 70
- ACCU-MEASURE® Critical Measurement Devices ..... 70
- RECEPTAL Solidifiers ..... 71
- AMSure® Suction Tubing with Connectors ..... 72
- Suction Yankauers and Kits ..... 73
- Suction Device Holster System ..... 74
- Accessories ..... 74
- Suction Catheters and Trays ..... 77
- Suction Catheter Kits ..... 78
- Reusable Hardware ..... 79
- Regulator Mounts ..... 82

## RECEPTAL® Outer Shells

RECEPTAL Outer Shells are designed with durability and efficiency in mind. The reusable canister is to be used with the same size RECEPTAL Disposable Liner.

- Easy to read volume marks
- Large patient label strip for convenient documentation
- Non-sterile, DEHP-free, PVC-free, not manufactured with natural rubber latex



Color-coded by size for easy selection.


REF	DESCRIPTION	QTY/UNIT
43449-05 ©	1000 mL Outer Shell	10/Case
43423-01 ©	1500 mL Outer Shell	10/Case
43445-01 ©	2000 mL Outer Shell	10/Case
43503-01 ©	3000 mL Outer Shell	8/Case



The canister ring seal firmly secures the liner lid to the canister, maximizing vacuum power to provide the best clinical outcomes.



Large, easy to read volume marks in 25 mL increments help with patient documentation.



Canister "T" designates vacuum source connection to simplify canister set up.

## RECEPTAL Liners

RECEPTAL Disposable Liners are designed with permanently attached lids and a multifunctional lid tube which help decrease the risk of fluid exposure. RECEPTAL Liners are to be used with the same size RECEPTAL Outer Shells.

- Anti-foaming emulsion helps prevent froth formation and aerosolization, reducing the risk of contaminant exposure
- 100% leak tested seamless liners minimize risk of contamination
- Non-sterile, BPA-free, not manufactured with natural rubber latex



SIZE (mL)	REF	DESCRIPTION	QTY/UNIT
BULK PACKAGED			
1000	43010-05	Liner with VAC-GARD® and Pour Spout	50/Case
	43010S ©	Liner with VAC-GARD and 30g Solidifier	50/Case
BULK PACKAGED			
1500	43025-01	Liner with Pour Spout	50/Case
	43044-01 ©	Liner with VAC-GARD and Pour Spout	50/Case
	43015S ©	Liner with VAC-GARD and 45g Solidifier	50/Case
	INDIVIDUALLY WRAPPED		
	0G903-01 ©	Liner with Pour Spout	50/Case
	0G904-01 ©	Liner with VAC-GARD	50/Case
2000	0L654-01 ©	Liner with VAC-GARD and Pour Spout	50/Case
	BULK PACKAGED		
	43027-01	Liner with Pour Spout,	50/Case
	43041-01	Liner with Ortho Port	50/Case
	43042-01	Liner with VAC-GARD	50/Case
	43044-05 ©	Liner with VAC-GARD and Pour Spout	50/Case
	43020S ©	Liner with VAC-GARD and 60g Solidifier	50/Case
	INDIVIDUALLY WRAPPED		
	0G901-01 ©	Liner	50/Case
	0G902-01 ©	Liner with VAC-GARD	50/Case
	0P508-01 ©	Liner with VAC-GARD and Pour Spout	50/Case
3000	0G905-01 ©	Liner with Ortho Port	50/Case
	0P507-01 ©	Liner with Pour Spout	50/Case
	BULK PACKAGED		
	0L656-01*	Liner with Dual Ports (Ortho + Standard) and Pour Spout	50/Case
0L655-01*	Liner with VAC-GARD, Dual Ports (Ortho + Standard) and Pour Spout	50/Case	
43030S* ©	Liner with VAC-GARD, Dual Ports (Ortho + Standard), Pour Spout and 90g Solidifier	50/Case	

\* Dual port allows for dual suctioning during procedures.  
Standard Port = 5/16" (7.9mm) OD, 1/4" (6.4mm) ID; Ortho Port = 3/8" (9.6mm) OD, 5/16" (7.8mm) ID

## RECEPTAL® Twist Liners

The RECEPTAL Twist liner is a patented all-in-one waste fluid collection and solidifier solution available in 1.5 and 2 L sizes, in both VAC-GARD® and non-VAC-GARD options.

After a procedure, simply twist the integrated and preloaded Solidifier Dispenser Cap at the top of the lid clockwise 180° to easily and safely release solidifier into the liner. The Twist eliminates the need to open ports, creating a closed system where waste fluid is safe, sealed, solid and ready for disposal.

- Premeasured, preloaded and contact-free Solidifier Dispenser Cap
- Integrated solidifier eliminates the need for open ports
- All-in-one solidifier and liner solution reduces costs associated with extra inventory
- 100% leak tested seamless liners minimize risk of contamination
- Non-sterile, BPA-free, not manufactured with natural rubber latex



43121S



43120S shown above in canister 43445-01 (canister not included)

Premeasured, preloaded and contact-free solidifier bottle

Anti-foaming emulsion helps prevent froth formation and aerosolization, reducing the risk of contaminant exposure

Designed with permanently attached lids and a multi-functional lid tube which helps decrease the risk of fluid exposure

VAC-GARD® option, a non-mechanical shut-off valve that supports best practice by eliminating 99.9% of pseudomonas and preventing backflow of aerosols and fluids into the vacuum system

REF	DESCRIPTION	QTY/UNIT
43110S	1000 mL Liner with VAC-GARD and Pre-attached Solidifier Dispenser Cap	40/Case
43115S ©	1500 mL Liner with VAC-GARD and Pre-attached Solidifier Dispenser Cap	40/Case
43116S ©	1500 mL Liner with Pre-attached Solidifier Dispenser Cap	40/Case
43120S ©	2000 mL Liner with VAC-GARD and Pre-attached Solidifier Dispenser Cap	40/Case
43121S ©	2000 mL Liner with Pre-attached Solidifier Dispenser Cap	40/Case



# iRECEPTAL® Surgical Waste Fluid Management System

The iReceptal Digital Surgical Suction System is an all-in-one, closed suction, container, and disposal system. This advanced, easy-to-use surgical suction solution is completely self-contained to provide for independent collection and disposal of waste fluids, along with a docking station for automatic and closed emptying of the rover.

**iReceptal is currently not cleared by the US FDA.**

- Constantly fully closed system, from collection to disposal
- Staff is never exposed to potentially infectious fluid waste
- Helps optimize surgical turnaround times
- Manifold provides multiple suction lines at the same time with backflow protection
- Separates fluid from tissue and bone fragments
- Waste fluid discharges directly into facility's drain system
- High volume capacity of 20L
- Volume gradation indications
- Easy to maneuver with large wheels
- Lowers disposal costs by eliminating full suction canisters in the waste stream
- Automatically cleans when docked in under 5 minutes

The system is comprised of:


- Mobile Rover
- Docking Station
- Disposable Suction Manifold
- Single-Use Smoke Filter
- Cleaning Detergent



Mobile Rover



Docking Station



### Manifold Storage

- 4 ports for suction tubing with backflow control



### Intelligent Control Panel

- Fluid read-out (Version 1.5 only)

### Completely Closed System

- Avoid direct touch of waste fluid

### Electromagnetic Couplers

- Allow for easy docking



### Direct to Drain Systems

- Hands-free disposal and self-cleaning

### 20L Volume

- Suitable for large amounts of fluid




REF	DESCRIPTION	QTY/UNIT
iRR100-01 ©	iRECEPTAL 1.0 Rover	1/Each
iRR150-01 ©	iRECEPTAL 1.5 Rover with 7cm (2.7") LCD Screen	1/Each
iRD102-01 ©	iRECEPTAL Docking Station (220VAC)	1/Each
iRD101-01 ©	iRECEPTAL Docking Station (110VAC)	1/Each
iRM100-01 ©	4-Port Manifold (disposable)	1/Each
iRC001-01 ©	Detergent for Docking Station	1/Each
iRC002-01 ©	Detergent for Rover	1/Each

## EZE-VAC® Suction Canister System

EZE-VAC suction canisters are a single-use disposable, designed for convenience by removing the need for a secondary waste liner. The canisters include VAC-GARD, a non-mechanical shut-off valve that supports best practice by eliminating 99.9% of pseudomonas and preventing backflow of aerosols and fluids into the vacuum system.

- Easy-to-read volume marks
- Reflux protector on patient port reduces contamination exposure risk
- Non-sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex



REF	DESCRIPTION	QTY/UNIT
43203-01	1200 mL Canister, Lid and 2 Elbows	32/Case
43206-01	2000 mL Canister, Lid and 2 Elbows	24/Case
43213-01	3300 mL Canister, Lid and 2 Elbows	12/Case

## EZE-VAC® 2.0 Suction Canister System

The new EZE-VAC 2.0 suction canister system is very economical and compatible with most hardware in the market. This canister system also includes VAC-GARD® non-mechanical shut-off valve and an anti-foaming emulsion that prevents froth formation and aerosolization, reducing the risk of contaminant exposure while maximizing fluid capacity.



- Easy-to-read volume marks
- VAC-GARD eliminates 99.9% of pseudomonas and protects the vacuum system from fluids
- Compatible with most hardware in the market
- Non-sterile, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	QTY/UNIT
43210-01	800 mL EZE-VAC 2.0 Canister and Lid	32/Case
43202-01	1200 mL EZE-VAC 2.0 Canister and Lid	32/Case
43202-820821	1200 mL EZE-VAC 2.0 Canister Procedure Kit with Non-sterile Suction Tube 3/16 in. ID x 6 ft. (5 mm x 1.8 m) + 3/16 in ID x 18 in. (5 mm x 46 cm)	30/Case

## ACCU-MEASURE® Critical Measurement Devices

ACCU-MEASURE critical measurement canisters enable the user to make contact-free critical measurements of smaller quantities of aspirate (i.e. during pediatric procedures).

The RECEPTAL® ACCU-MEASURE is a disposable, tapered 250 mL critical measurement canister calibrated to scale. Graduations in 5 mL increments from 0 to 50 mL, 10 mL increments from 50 mL to 150 mL and 20 mL increments from 150 mL to 250 mL.

- Non-sterile, BPA-free, not manufactured with natural rubber latex



REF	DESCRIPTION	FEATURES	QTY/UNIT
43306-01	250 mL Critical Measurement Canister	Tandem Tube	12/Case
43305-10	Canister Attachment Bracket for 250 mL Canister	Compatible with 1000 mL, 1500 mL and 2000 mL RECEPTAL Canisters	10/Case



## RECEPTAL® Solidifier Bottles and Dissolvable Pouches

RECEPTAL Solidifier Bottles and Dissolvable Pouches effectively solidify potentially serious infectious waste within the suction canister. It penetrates rapidly throughout the canister/liner to solidify blood and body fluids and gels. It retains solidity to help prevent spills and leaks for safe and easy bagging and transport to disposal.

- No shaking required
- Multiple sizes to match any need
- Dissolvable pouch option eliminates plastic bottle waste
- Chlorine-free for non-toxic incineration
- Not made with natural rubber latex



4315SB

4320SB

### RECEPTAL Solidifier Bottles

REF	DESCRIPTION	QTY/UNIT
<b>4305SB</b>	500 Solidifier Bottle, solidifies 500 mL	96/Case
<b>4312SB</b>	1200 Solidifier Bottle, solidifies 1200 mL	96/Case
<b>4315SB</b>	1500 Solidifier Bottle, solidifies 1500 mL	96/Case
<b>4320SB</b>	2000 Solidifier Bottle, solidifies 2000 mL	96/Case
<b>4330SB</b>	3000 Solidifier Bottle, solidifies 3000 mL	96/Case
<b>43160SB</b>	16000 Solidifier Bottle, solidifies up to 16,000 mL	20/Case



4320SP

### RECEPTAL Solidifier Pouches

REF	DESCRIPTION	QTY/UNIT
<b>4305SP</b>	500 Solidifier Pouch, solidifies 500 mL	96/Case
<b>4312SP</b>	1200 Solidifier Pouch, solidifies 1200 mL	96/Case
<b>4315SP</b>	1500 Solidifier Pouch, solidifies 1500 mL	96/Case
<b>4320SP</b>	2000 Solidifier Pouch, solidifies 2000 mL	96/Case
<b>4330SP</b>	3000 Solidifier Pouch, solidifies 3000 mL	96/Case

## AMSure® Suction Tubing with Connectors

AMSure Suction Tubing, with molded female connectors and male adapters, provide secure connections to support uninterrupted procedures. Flexible, low-coil memory allows for easy draping.

- Sterile and non-sterile options, BPA-free, not manufactured with natural rubber latex



AS825

DIAMETER (ID)	LENGTH	STERILE	NON-STERILE
3/16 in. (5 mm) Green Connectors	18 in. (46 cm)	<b>AS820</b> (50/Case)	<b>AS820B</b> (100/Case)
	6 ft. (1.8 m)	<b>AS821</b> (50/Case)	<b>AS821B</b> (50/Case)
	10 ft. (3.0 m)	<b>AS822</b> (50/Case)	<b>AS822B</b> (50/Case)
	12 ft. (3.7 m)	<b>AS823</b> (20/Case)	<b>AS823B</b> (50/Case)
	20 ft. (6.1 m)	<b>AS824</b> (20/Case)	<b>AS824B</b> (50/Case)
1/4 in. (6 mm) Light Blue Connectors	6 ft. (1.8 m)	<b>AS825</b> (50/Case)	<b>AS825B</b> (50/Case)
	10 ft. (3.0 m)	<b>AS826</b> (50/Case)	<b>AS826B</b> (50/Case)
	10 ft. (6 mm x 3.0 m) with Pinch Clamps	<b>AS826C</b> (50/Case)	—
	12 ft. (3.7 m)	<b>AS827</b> (20/Case)	<b>AS827B</b> (50/Case)
	20 ft. (6.1 m)	<b>AS828</b> (20/Case)	<b>AS828B</b> (50/Case)
	100 ft. (30.0 m) with Connectors every 6 ft. (1.8 m)	—	<b>AS829B-01</b> (1Roll/Case)
9/32 in. (7 mm) Dark Blue Connectors	6 ft. (1.8 m)	<b>AS825M</b> (50/Case)	<b>AS825BM</b> (50/Case)
	10 ft. (3.0 m)	<b>AS826M</b> (50/Case)	<b>AS826BM</b> (50/Case)
	12 ft. (3.7 m)	<b>AS827M</b> (20/Case)	<b>AS827BM</b> (50/Case)
	20 ft. (6.1 m)	<b>AS828M</b> (20/Case)	<b>AS828BM</b> (50/Case)

## Suction Yankauers and Kits

AMSure® Suction Yankauers are designed with a slip-resistant, lightweight handle that helps provide optimal control and maneuverability.

- Translucent pathway allows for fluid visualization
- Smooth tips are designed with side eyes to minimize the chance of tissue grab and damage
- Sterile, disposable, DEHP-free, BPA-free, not manufactured with natural rubber latex



### Yankauers

DESCRIPTION			STERILE	NON-STERILE
Rigid	Regular Beveled Tip	Non-Vented	<b>AS830</b> © 50/Case	<b>AS830B</b> 500/Case
Rigid	Regular Beveled Tip	Vented	<b>AS831</b> © 50/Case	<b>AS831B</b> 500/Case
Rigid	Bulbous Tip	Non-Vented	<b>AS832</b> © 50/Case	<b>AS832B</b> 500/Case
Rigid	Bulbous Tip	Vented	<b>AS833</b> © 50/Case	<b>AS833B</b> 500/Case
Flexible	Regular Beveled Tip	Non-Vented	<b>AS834</b> © 50/Case	<b>AS834B</b> 500/Case
Flexible	Regular Beveled Tip	Vented	<b>AS835</b> © 50/Case	<b>AS835B</b> 500/Case
Flexible	Bulbous Tip	Non-Vented	<b>AS836</b> © 50/Case	<b>AS836B</b> 500/Case
Flexible	Bulbous Tip	Vented	<b>AS837</b> © 50/Case	<b>AS837B</b> 500/Case

### Yankauer & Tubing Kits

REF	DESCRIPTION	QTY/UNIT
<b>AS825P</b> ©	Rigid Yankauer, Bulbous Tip, Non-Vented, pre-attached to 1/4 in. (6 mm) × 6 ft. (1.8 m) Suction Tubing, Sterile	50/Case
<b>AS827P*</b> ©	Rigid Yankauer, Non-Vented with Bulbous Tip, pre-attached to 1/4 in. (6 mm) × 12 ft. (3.7 m) Suction Tubing, Sterile	50/Case

\*Made to order

## Suction Device Holster System

SteriCadi™ is a sterile disposable holder to store suction devices. The intuitive design is easy to set up and use to help support better infection control outcomes.

- Suction devices can be isolated to help reduce cross-contamination and improve safety measures
- Sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	QTY/UNIT
<b>SIH0001</b>	SteriCadi Suction Device Holder with Clip	30/Case



## Accessories

### RECEPTAL® Specimen Traps

RECEPTAL Specimen Traps safely collect and retain specimens for analysis. The traps are compatible with RECEPTAL and EZE-VAC® lids with a pour spout.

- Pour spout cover provides a tight seal to ensure proper suction
- Specimen traps with stiffener offer simplified insertion into the collection port
- Non-sterile, disposable, DEHP-free, BPA-free, not manufactured with natural rubber latex



REF	DESCRIPTION	FEATURES		QTY/UNIT
		Suction Port Size	Stiffener	
<b>43303-01</b> ©	4 in. (10.16 cm), Pour Spout Cover and Tethered Cap	3/8 in. (0.95 cm)		24/Case
<b>43303-02</b> ©	4 in. (10.16 cm), Pour Spout Cover and Tethered Cap	3/8 in. (0.95 cm)	•	24/Case
<b>43303-05</b> ©	8 in. (20.32 cm), Pour Spout Cover and Tethered Cap	1/2 in. (1.27 cm)	•	24/Case

### External VAC-GARD® Shut-Off Valve

VAC-GARD external non-mechanical shut-off valve helps support best practice by eliminating 99.9% of pseudomonas and preventing backflow of aerosols and fluids into the vacuum system.

- Non-sterile, disposable, DEHP-free, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	QTY/UNIT
<b>43309-01</b>	External Shut-Off Valve, 1/4 in. (0.64 cm) Male × Female Elbow	10/Case





## Accessories



### RECEPTAL® Accu-Measure® Device Bracket

This bracket is used to secure 250 mL canister (43306-01) to a RECEPTAL Reusable Canister. Compatible with 1000 mL, 1500 mL and 2000 mL RECEPTAL Canisters.

REF	DESCRIPTION	QTY/UNIT
43305-10	Canister Attachment Bracket for 250 mL Canister	10/Case

### Suction Tubing

For permanent Outer Shell-to-vacuum source, or Outer Shell-to-Outer Shell connections.



REF	DESCRIPTION	QTY/UNIT
43455-01	100 ft. (30.48 m) PVC Non-Sterile Suction Tubing 9/32 in. (7 mm) I.D.	1/Case
43457-25	24 in. (0.61 m) PVC Non-Sterile Suction Tubing 9/32 in. (7 mm) I.D.	25/Case

### EZE-VAC® Tandem Tubing

EZE-VAC Tandem disposable tubing allows two or more EZE-VAC canisters to work together to meet procedural requirements.


- Non-sterile, BPA-free, not manufactured with natural rubber latex



REF	DESCRIPTION	QTY/UNIT
43317-25	18 in. (0.45 m) Tandem Tubing	25/Case

### RECEPTAL Outer Shell Connection Tubing


RECEPTAL Outer Shell connection disposable tubing allows to easily connect the vacuum "T" of the RECEPTAL canisters between them for easy tandem connection.



REF	DESCRIPTION	QTY/UNIT
CTM50B ©	1/4 in. x 20 in. (6.35 mm x 500 mm), Female Connector	100/Case
CTM20B ©	1/4 in. x 8 in. (6.35 mm x 200 mm), Female Connector	300/Case

### Tubing Manifold

For connection of 2, 3 or 4 canisters.



REF	DESCRIPTION	QTY/UNIT
TMF001	Tubing Manifold	10/Case

## RECEPTAL® Elbow Adapters and Taper

Use with Patient Port.

REF	DESCRIPTION	QTY/UNIT
43318-01 ©	1/4 in. (0.64 cm) Male x Female Elbow (use with 3/8 in. (0.95 cm) Orthopedic Port)	100/Case
43316-01 ©	1/4 in. (0.64 cm) Male x Female Elbow	100/Case
43494-25	1/4 in. (0.64 cm) Male x Female Elbow	25/Case
43473-25	Double 1/4 in. (0.64 cm) Male Taper (use with 3/8 in. (0.95 cm) Orthopedic Port)	25/Case



### RECEPTAL Outer Shell Parts

Attach to Canister for vacuum suction.

REF	DESCRIPTION	QTY/UNIT
44300-05	T, Double Male 1/4 in. (0.64 cm) with O-Ring and Screw	5/Case
43448-25	Taper, Male 1/4 in. (0.64 cm) x Female (use with Canister Vacuum "T")	25/Case



### 5-in-1 Straight Connector

AMSure® 5-in-1 Connector is a durable connector used to connect same size tubing or use as a reducer and fits five (3/16 in., 1/4 in., 5/16 in., 3/8 in., 7/17 in.) sizes of tubing.

- DEHP-free, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	QTY/UNIT
AS807 ©	5-in-1 Connector, Sterile	100/Case
AS807B	5-in-1 Connector, Non-sterile	500/Case



### 6-in-1 Universal Y Connector

Lightweight, durable polypropylene plastic connectors adapt for almost any use allowing for multiple pieces of tubing to be joined together. Tapered legs fit six sizes of tubing, from 6 mm to 15 mm (1/4 in. to 9/16 in.).

- DEHP-free, BPA-free, not manufactured with natural rubber latex

REF	DESCRIPTION	QTY/UNIT
UYC001S	6-in-1 Universal Y Connector, Sterile	100/Case
UYC001	6-in-1 Universal Y Connector, Non-sterile	500/Case



## Suction Catheters and Trays

### Suction Catheters

AMSure® Suction Catheters improve patient care by safely removing respiratory secretions from the upper airway. The suction catheters are available in both straight and coiled packs in various French sizes.

- **Opposing staggered eyelets reduce the potential for blocked suction**
- **Whistle cap port enhances controlled suction**
- **Sterile, BPA-free, not manufactured with natural rubber latex**

AS361



REF	DESCRIPTION	FEATURES	QTY/UNIT
<b>AS361</b> ©	6 Fr, 17 ½ in. (44.5 cm)	Straight Pack	50/Case
<b>AS361C</b> ©	6 Fr, 17 ½ in. (44.5 cm)	Coiled Pack	50/Case
<b>AS362</b> ©	8 Fr, 17 ½ in. (44.5 cm)	Straight Pack	50/Case
<b>AS362C</b> ©	8 Fr, 17 ½ in. (44.5 cm)	Coiled Pack	50/Case
<b>AS363</b> ©	10 Fr, 23 in. (58.4 cm)	Straight Pack	50/Case
<b>AS363C</b> ©	10 Fr, 23 in. (58.4 cm)	Coiled Pack	50/Case
<b>AS364</b> ©	12 Fr, 23 in. (58.4 cm)	Straight Pack	50/Case
<b>AS364C</b> ©	12 Fr, 23 in. (58.4 cm)	Coiled Pack	50/Case
<b>AS365</b> ©	14 Fr, 23 in. (58.4 cm)	Straight Pack	50/Case
<b>AS365C</b> ©	14 Fr, 23 in. (58.4 cm)	Coiled Pack	50/Case
<b>AS366</b> ©	16 Fr, 23 in. (58.4 cm)	Straight Pack	50/Case
<b>AS366C</b> ©	16 Fr, 23 in. (58.4 cm)	Coiled Pack	50/Case
<b>AS367</b> ©	18 Fr, 23 in. (58.4 cm)	Straight Pack	50/Case
<b>AS367C</b> ©	18 Fr, 23 in. (58.4 cm)	Coiled Pack	50/Case

### Mini Suction Catheter Trays

The AMSure Mini Suction Catheter Trays' space-saving design offers everything needed for effective procedures. The mini tray serves as a basin to hold sterile water or saline.

- **Opposing staggered eyelets reduce the potential for blocked suction**
- **DeLee Tip's smooth and rounded tip minimizes trauma for pediatric patients**
- **Whistle tip enhances controlled suction**
- **Sterile, BPA-free, not manufactured with natural rubber latex**



SCT10

REF	DESCRIPTION	CATHETER SIZE	FEATURES	QTY/UNIT
<b>SCT08</b>	Tray, Vinyl Gloves and 17 ½ in. Suction Catheter	8 Fr	DeLee Tip	100/Case
<b>SCT10</b>	Tray, Vinyl Gloves and 23 in. Suction Catheter	10 Fr	DeLee Tip	100/Case
<b>SCT12</b>	Tray, Vinyl Gloves and 23 in. Suction Catheter	12 Fr	Whistle Tip	100/Case
<b>SCT14</b>	Tray, Vinyl Gloves and 23 in. Suction Catheter	14 Fr	Whistle Tip	100/Case
<b>SCT16</b>	Tray, Vinyl Gloves and 23 in. Suction Catheter	16 Fr	Whistle Tip	100/Case

## Suction Catheter Kits

AMSure® Suction Catheter Kits provide hassle-free and convenient extraction of mucus or saliva from the upper respiratory airway.

- **Uniquely designed two-eyed whistle tip offers safe and effective suctioning**
- **Sterile, BPA-free, not manufactured with natural rubber latex**



AS371






REF	DESCRIPTION	CATHETER SIZE	FEATURES	QTY/UNIT
<b>AS371</b> ©	17 ½ in. (44.5 cm) Whistle Tip Catheter	6 Fr	1 Vinyl Glove	50/Case
<b>AS381</b> ©	17 ½ in. (44.5 cm) Whistle Tip Catheter	6 Fr	2 Vinyl Gloves	50/Case
<b>AS372</b> ©	17 ½ in. (44.5 cm) Whistle Tip Catheter	8 Fr	1 Vinyl Glove	50/Case
<b>AS382</b> ©	17 ½ in. (44.5 cm) Whistle Tip Catheter	8 Fr	2 Vinyl Gloves	50/Case
<b>AS373</b> ©	23 in. (58.4 cm) Whistle Tip Catheter	10 Fr	1 Vinyl Glove	50/Case
<b>AS383</b> ©	23 in. (58.4 cm) Whistle Tip Catheter	10 Fr	2 Vinyl Gloves	50/Case
<b>AS374</b> ©	23 in. (58.4 cm) Whistle Tip Catheter	12 Fr	1 Vinyl Glove	50/Case
<b>AS384</b> ©	23 in. (58.4 cm) Whistle Tip Catheter	12 Fr	2 Vinyl Gloves	50/Case
<b>AS375</b> ©	23 in. (58.4 cm) Whistle Tip Catheter	14 Fr	1 Vinyl Glove	50/Case
<b>AS385</b> ©	23 in. (58.4 cm) Whistle Tip Catheter	14 Fr	2 Vinyl Gloves	50/Case
<b>AS376</b> ©	23 in. (58.4 cm) Whistle Tip Catheter	16 Fr	1 Vinyl Glove	50/Case
<b>AS386</b> ©	23 in. (58.4 cm) Whistle Tip Catheter	16 Fr	2 Vinyl Gloves	50/Case
<b>AS387</b> ©	23 in. (58.4 cm) Whistle Tip Catheter	18 Fr	2 Vinyl Gloves	50/Case





## Reusable Hardware

### Canister Brackets

REF	DESCRIPTION	QTY/UNIT
 <b>43444-05</b>	Wall channel bracket for mounting canister, stainless steel	5/Case
 <b>43444P</b>	Wall channel bracket for mounting canister, polycarbonate	50/Case
 <b>44303-10</b>	Bracket Clamp Adapter for adapting Wall Bracket (43444-05) to fit stands, carts or vertical poles, stainless steel	10/Case
 <b>44304-05</b>	Canister Support Bracket to attach Canister directly to a vacuum regulator, use with 1000 mL, 1200 mL and 1500 mL canisters	5/Case
 <b>RCB001</b>	Metal Rail Bracket with Universal Holder Clamp	30/Case

## Regulator Mounts

REF	DESCRIPTION	QTY/UNIT
<b>43408-05</b>	Use to attach regulator alone or to other supports. 1/8 in. (0.32 cm) male pipe thread × 1/4 in. (0.64 cm) male taper with 1-1/8 in. (2.86 cm) threaded shank and 2 locknuts	5/Case
<b>43422-05</b>	Use when regulator is wall mounted to vacuum and a simple taper is needed to attach tubing to regulator. 1/8 in. (0.32 cm) male pipe thread × 1/4 in. (0.64 cm) male taper	5/Case
<b>43401-05</b>	Use to attach regulator alone or to other supports. 1/4 in. (0.64 cm) taper with 1-1/4 in. (3.18 cm) shank	5/Case
<b>43400-05</b>	Use when regulator is wall mounted to vacuum and a simple taper is needed to attach tubing to regulator. 1/4 in. (0.64 cm), male taper	5/Case
<b>OFV001</b>	Metal on-off valve, 1/8 in. NPT	25/Case

Note: Amsino does not sell Regulators





# Index

REF	PAGE	REF	PAGE	REF	PAGE	REF	PAGE	REF	PAGE	REF	PAGE	REF	PAGE	REF	PAGE
000606	16	43316-01	76	AS011	4, 39	AS376	49, 78	AS75035	46	AS85012	33	ASR0035E	48, 54	MR4207	16
003201	16	43317-25	75	AS015	39	AS381	49, 78	AS75040	44	AS85012K	33	ASR0039	56	MR4307	17
003202	16	43318-01	76	AS016	39	AS382	49, 78	AS75050	44	AS85014	33	ASR0039E	56	MR4B07	17
0G901-01	65	43400-05	82	AS032L	60	AS383	49, 78	AS75060	45	AS85014K	33	ASR0055	48, 54	MR4P07	18
0G902-01	65	43401-05	82	AS032M	60	AS384	49, 78	AS75080	43	AS85016	33	ASR0055E	48, 54	MR4R10	16
0G903-01	65	43408-05	82	AS032S	60	AS385	49, 78	AS75085	43	AS85016K	33	ASR0059	56	MR4T07	16
0G904-01	65	43422-05	82	AS0352	48, 54	AS386	49, 78	AS75180	44	AS860	47	ASR0059E	56	MR4TF04	17
0G905-01	65	43423-01	64	AS0352E	48, 54	AS387	49, 78	AS75181	43	AS860608	35	ASR0159	56	MR4Y32	16
0L654-01	65	43424-10	80	AS0552	48, 54	AS398	56	AS75190	44	AS860610	35	ASR0159E	56	MRA102-84	9
0L655-01	65	43424-V	80	AS0552E	48, 54	AS398E	56	AS75199	44	AS860612	35	AY0200	18	MRA10E-104	8
0L656-01	65	43425-05	80	AS1000	48, 54	AS399	56	AS75280	43	AS860614	35	CTM20B	75	MRA10E-95	7
0P507-01	65	43426-05	80	AS1065	48, 54	AS399E	56	AS75282	43	AS860616	35	CTM50B	75	MRA15-100	10
0P508-01	65	43426P	80	AS1065E	48, 54	AS409	37	AS75285	43	AS861	47	D38301	14	MRA15B-105	11
1010001B	9	43428-05	80	AS1066	48, 54	AS40008	29	AS76007	46	AS861006	35	D38302	14	MRA60B-92	13
107201	7	43428-V	80	AS115	4, 39	AS40010	29	AS78010	45	AS861008	35	D38305	14	MRS60B-106	14
108301	7	43429-01	80	AS116	4, 39	AS41012S	30	AS78015	45	AS861010	35	E-0500M	3	MRS60B-117	14
108302	9	43431-01	79	AS121	5, 38	AS41014	29	AS78020	45	AS861608	35	E-1004	3	MRSA100-84	7
108303	7	43432-01	79	AS126	5, 38	AS41014H	28	AS78025	45	AS861610	35	E-1200	3	MRSA103-94	9
108304	7	43436-05	80	AS127	5, 38	AS41014S	30	AS807	62, 76	AS861610C	34	E-1203	3	MRSA103-94	9
108305	8	43438-10	80	AS130	5, 38	AS41016	29	AS807B	62, 76	AS861612	35	ENF1004K	3	MRSA10E-107	8
108306	7	43438-V	80	AS136	5, 38	AS41016H	28	AS820	51, 58, 72	AS861612C	34	ENF1200K	3	MRSA15-100	10
108308	8	43439-01	80	AS1565	48, 54	AS41016S	30	AS820B	51, 58, 72	AS861614	35	ENF1203K	3	MRSA15-91	10
109602	7	43444-05	81	AS1565E	48, 54	AS41018	29	AS821	51, 58, 72	AS861614C	34	ENF500K	3	MRSA15E-120	10
1510501	10	43444P	81	AS1566	48, 54	AS41018H	28	AS821B	51, 58, 72	AS861616	35	ENK127	5	MRSA15E-88	10
158305	10	43445-01	64	AS1765	48, 54	AS41018S	30	AS822	51, 58, 72	AS861616C	34	ENK136	5	MRSA15E-96	11
158306	10	43447-01	79	AS1765E	48, 54	AS41020	29	AS822B	51, 58, 72	AS861618	35	ENS015	4	MRSA60-100	13
158307	10	43448-25	76	AS302	26	AS41020S	30	AS823	51, 58, 72	AS861618C	34	ENS016	4	MRSA60-84	13
158309	10	43449-05	64	AS306N	37	AS41022	29	AS823B	51, 58, 72	AS870	31	ENS115	4	MRSA602-100	13
158310	9	43455-01	75	AS307N	37	AS41022S	30	AS824	51, 58, 72	AS87114	31	ENS115NS	4	MRSA60E-88	13
159602	10	43457-25	75	AS308N	37	AS41024	29	AS824B	51, 58, 72	AS87118	31	ENS116	4	MRSB-102	15
43010-05	65	43467-01	80	AS309N	37	AS41024S	30	AS825	51, 58, 72	AS880	25	FP092	19	MRSB-105	15
43010S	65	43469-12	80	AS3100	11	AS41026	29	AS825B	51, 58, 72	AS880K	25	FP092Y1	19	MRSB-94	15
43015S	65	43469-V	76	AS3101	11	AS41026S	30	AS825BM	43, 50, 60	AS890	25	FP100B12	20	MRSFC-19	17
43020S	65	43473-25	76	AS3102	11	AS41028	29	AS825M	51, 58, 72	AS890K	25	FP100BY2	19	OFV001	82
43025-01	65	43474-01	79	AS3110A	11	AS41028S	30	AS825P	52, 59, 73	AS89214	32	FP100FY1	19	OM001	42
43027-01	65	43475-01	79	AS3110B	11	AS41030	29	AS826	51, 58, 72	AS89214R	32	FP106F12	20	PE90B12T	2
43030S	65	43476-P2	79	AS312	26	AS41032	29	AS826B	51, 58, 72	AS89214RK	32	FP106FY2	19	PE90B5T	2
43041-01	65	43476-P3	79	AS322	26	AS42012	29	AS826BM	51, 58, 72	AS89216RK	32	FP109B13	20	PE90ST	2
43042-01	65	43494-25	76	AS32200	26	AS42012S	30	AS826C	51, 58, 72	AS89300	23	FP109BY3	19	PG001	2
43044-01	65	43503-01	64	AS326	26	AS42014	29	AS826M	51, 58, 72	AS89314	22	FP169BY3	19	RCB001	81
43044-05	65	44300-05	76	AS328	27	AS42014H	28	AS827	51, 58, 72	AS89314H	22	HME001	47	S-1000-SI	3
43055B	71	44303-10	81	AS32600	26	AS42014S	30	AS827B	51, 58, 72	AS89314S	22	HME002	47	SCT08	50, 77
43055P	71	44304-05	81	AS327100	27	AS42016H	28	AS827BM	51, 58, 72	AS89316	22	iRC001-01	67	SCT10	50, 77
43110S	66	607201	13	AS330	61	AS42016	29	AS827M	51, 58, 72	AS89316H	22	iRC002-01	67	SCT12	50, 77
43115S	66	608304	13	AS333	61	AS42016S	30	AS827P	52, 59, 73	AS89316S	22	iRD101-01	67	SCT14	50, 77
43116S	66	608305	13	AS340	40, 60	AS42018	29	AS828	51, 58, 72	AS89318	22	iRD102-01	67	SCT16	50, 77
43120S	66	608306	13	AS341	40, 60	AS42018H	28	AS828B	51, 58, 72	AS89318H	22	iRM100-01	67	SIH0001	59, 74
43121S	66	AA3101	9	AS343	40, 60	AS42018S	30	AS828BM	51, 58, 72	AS89318S	22	iRR100-01	67	TMF001	75
43125B	71	ABL400	36	AS345	40, 60	AS42020	29	AS828M	51, 58, 72	AS89330	23	iRR150-01	67	UCC100	60
43125P	71	ABL402	36	AS3500	15	AS42020S	30	AS829B-01	51, 58, 72	AS89800	25	IVF1003	55	UYC001	62, 76
43155B	71	ABL404	36	AS3515	15	AS42022	29	AS830	52, 59, 73	AS89814	24	IVF1005	55	UYC001S	62, 76
43155P	71	AC3110	7	AS361	49, 77	AS42024	29	AS830B	52, 59, 73	AS89814H	24	IVF1010	55	VS363RC124	50
43160SB	71	AE0106	17	AS361C	49, 77	AS42026	29	AS831	52, 59, 73	AS89814S	24	IVF1010T	55	VS363TC124	50
43202-01	70	AE0107	17	AS362	49, 77	AS42028	29	AS831B	52, 59, 73	AS89816	24	IVF1010TM	55	VS363TP124	50
43202-820821	70	AE1106	17	AS362C	49, 77	AS42030	29	AS832	52, 59, 73	AS89816H	24	IVF1203	55	VS363TTP124	50
43203-01	69	AE1112	17	AS363	49, 77	AS44008	35	AS832B	52, 59, 73	AS89816S	24	IVF1203SS	55	VS364RC124	50
43206-01	69	AE1138	17	AS363C	49, 77	AS44010	35	AS833	52, 59, 73	AS89818	24	IVF1205	55	VS364TC124	50
43205B	71	AE3107	16	AS364	49, 77	AS44012	35	AS833B	52, 59, 73	AS89818H	24	IVF1205SS	55	VS364TP124	50
43205P	71	AE3108	16	AS364C	49, 77	AS44014	35	AS834	52, 59, 73	AS89818S	24	IVF1210	55	VS364TTP124	50
43210-01	70	AFS102	8	AS365	49, 77	AS44016	35	AS834B	52, 59, 73	AS89830	25	IVF1210SS	55	VS365RC124	50
43213-01	69	AFS107	17	AS365C	49, 77	AS44018	35	AS835	52, 59, 73	AS960610	34	IVF2020	55	VS365TC124	50
43303-01	74	AFS208	8	AS366	49, 77	AS44020	35	AS835B	52, 59, 73	AS960612	34	IVF2020SS	55	VS365TP124	50
43303-02	74	AFS208B	8	AS366C	49, 77	AS44022	35	AS836	52, 59, 73	AS960614	34	LBE306N	37	VS365TTP124	50
43303-05	74	AS00501	61	AS367	49, 77	AS44024	35	AS836B	52, 59, 73	AS960616	34	LBE307N	37		
43305-10	70	AS00501S	61	AS367C	49, 77	AS44028	35	AS837	52, 59, 73	AS961010	34	LBE308N	37		
43305-10	75	AS00502	61	AS371	49, 78	AS730100	46	AS837B	52, 59, 73	AS961610	34	LBE309N	37		
43306-01	70	AS00502S	61	AS372	49, 78	AS730200	46	AS85008	33	AS961612	34	MB4V01	18		
43309-01	74	AS00503	61	AS373	49, 78	AS74010	42	AS85008K	33	AS961614	34	MBSS001	18		
43305B	71	AS00503S	61	AS374	49, 78	AS75010	44	AS85010	33	AS961616	34	MR4001	18		
43305P	71	AS010	4, 39	AS375	49, 78	AS75020	44	AS85010K	33	ASR0035	48, 54	MR4107	16		



Advancing Healthcare Worldwide<sup>®</sup>

**Amsino Medical Group Co Ltd.**

708 Corporate Center Drive  
Pomona, CA 91768  
USA

P: +1 800.632.6746

F: +1 909.626.3888

info@amsino.com

**[www.amsino.com](http://www.amsino.com)**

---

All products listed are legally marketed and sold within the US and, where required, cleared by the US FDA. Products listed may not be available or licensed in all geographies. Amsino's manufacturing facilities are certified to ISO 13485:2016 and MDSAP. Please contact our customer service department or our regional office for product availability in your country or region.

Items listed with © are CE marked and legally available within the European Union (EU) and other regions where CE marking may be required.

---

Product(s) described may not be licensed or available for sale in all countries. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. AMSINO®, ACCU-MEASURE®, AMSafe®, AMSafe-3®, AMShield®, AMSure®, AMSure HG®, EZE-VAC®, IRECEPTAL®, One Mask™, Puggle®, R-Polished®, Receptal®, Safgard™, SteriCadi®, Sure-Lok™, THE-POLE-SYRINGE® and VAC-GARD® are registered trademarks of Amsino Medical Group Co. Ltd. © Copyright 2023, Amsino Medical Group Co. Ltd. All rights reserved.