

Urological Care Product Catalog





Amsino Medical Group offers a comprehensive range of urological products, for acute care, long-term care, and home care. Whether short-term or for chronic care, Amsino's urological care products are designed for optimal patient comfort, safety, and infection prevention.

Table of Contents

Foley Catheter Trays & Collection	1
Pre-Connected Foley Catheter Urine Drainage Bag Trays	1
Add-a-Foley Catheter Urine Drainage Bag Trays	2
Pre-Connected Foley Catheter Urine Meter Kits	3
Add-a-Foley Catheter Urine Meter Trays	4
Foley Insertion Trays	5
Urine Drainage Bags	7
Urine Meter Drainage Bags	9
Foley Catheters	11
Hydrophilic-Coated Latex	11
Silicone-Coated Latex	12
100% Silicone	13
Intermittent Urinary Catheters	14
Urethral Insertion Trays and Kits	14
Closed System Urethral Drainage Trays	15
Self-Catheterization Kits with R-Polished® Eyes	16
PVC Intermittent Catheters with R-Polished® Eyes	17
PVC Intermittent Urethral Catheter	18
Red Rubber Latex Urethral Catheters	19
Irrigation Sets	20
Ancillary & Patient Care Products	21
Leg Bags	21
Leg Bag Extension Tubes	21
Pediatric Urine Collector	22
Irrigation Trays & Kits	23
Syringes	24
Specimen Containers	25
Index of Item Codes	26

Foley Catheter Trays & Collection



AS89316

AMSure® Pre-Connected Foley Catheter Urine Drainage Bag Trays

AMSure Pre-Connected Foley Catheter Urine Drainage Bag Trays improve simplicity and efficiency of care. The catheter is pre-attached to a low-profile drainage bag which helps to minimize contamination. **(HCPCS codes: See product chart below)**

- Tray includes all necessary components for aseptic catheter insertion
- Tray improves clinical efficacy through cost-effective inventory management
- Gloves and povidone iodine solution help facilitate aseptic technique and reduce the risk of infection
- Needle-free urinary sampling port reduces the risk of needle stick injuries
- Trays are available with 100% silicone, silicone-coated latex, or hydrophilic-coated latex Foley catheters to suit patient needs
- Sterile

Trays include: 2000 mL low-profile urinary drainage bag, pre-connected 16 in. Foley catheter with 5 cc balloon, waterproof and fenestrated drapes, Nitrile powder-free gloves, (5) Rayon balls, 10% Povidone Iodine solution, forceps, lubricating jelly, syringe prefilled with 10 cc sterile water, 2 oz. specimen container with a patient label, central supply room wrap and Tyvek header bag.

CATHETER TYPE	ITEM CODE	FOLEY SIZE	BALLOON SIZE	QTY/UNIT
100% Silicone (HCPCS Code A4315)	AS89314S	14 Fr.	5 cc	10/Case
	AS89316S	16 Fr.	5 cc	10/Case
	AS89318S	18 Fr.	5 cc	10/Case
Silicone Coated Latex (HCPCS Code A4314)	AS89314	14 Fr.	5 cc	10/Case
	AS89316	16 Fr.	5 cc	10/Case
	AS89318	18 Fr.	5 cc	10/Case
Hydrophilic Coated Latex (HCPCS Code A4314)	AS89314H	14 Fr.	5 cc	10/Case
	AS89316H	16 Fr.	5 cc	10/Case
	AS89318H	18 Fr.	5 cc	10/Case



AS89300



AMSure® Add-a-Foley Catheter Trays

AMSure Add-a-Foley Catheter Trays improve simplicity and efficiency of care by including the urine collection device and necessary supplies to safely insert a Foley catheter. **(HCPCS code A4354)**

- Needle-free sampling port on urinary bag reduces the risk of needle stick injuries
- Gloves and povidone iodine swab sticks help maintain sterility and minimized risk of infection
- Sterile, not manufactured with natural rubber latex

Trays include: 2000 mL low-profile urinary drainage bag, waterproof and fenestrated drapes, Nitrile powder-free gloves, (5) Rayon balls, 10% povidone iodine solution, forceps, lubricating jelly, 2 oz. specimen container with patient label, central supply room wrap and Tyvek header bag.

ITEM CODE	DESCRIPTION	PRE-FILLED SYRINGE	QTY/UNIT
AS89300	Add-a-Foley Tray 2000 mL Urine Drainage Bag	10 cc	10/Case
AS89330	Add-a-Foley Tray 2000 mL Urine Drainage Bag	30 cc	10/Case

**BZK
Swabstick
Kits**

COMING SOON!

AMSure® Pre-Connected Foley Catheter Urine Meter Kits

The AMSure pre-connected urine meter trays offer the benefit of a closed system along with the accuracy of a meter to measure exact urine output. These truly closed systems provide all the necessary supplies for an aseptic catheter insertion to promote the US Centers for Disease Control (CDC) recommendations to help prevent catheter associated urinary tract infections (CAUTIs).

(HCPCS codes: See product chart below)

Trays include: 400 ml Urine Meter Preconnected with 5 cc Foley Catheter, waterproof and fenestrated drapes, (2) Nitrile gloves, 10% Povidone Iodine solution, (5) Rayon balls, forceps, lubricating jelly, syringe prefilled with 10 cc sterile water, 2 oz. specimen container with a patient label, central supply room wrap, Tyvek header bag, Sterile, not manufactured with rubber latex.



CATHETER TYPE	ITEM CODE	FOLEY SIZE	BALLOON SIZE	QTY/UNIT
100% Silicone (HCPCS Code A4315)	AS89814S	14 Fr.	5 cc	10/Case
	AS89816S	16 Fr.	5 cc	10/Case
	AS89818S	18 Fr.	5 cc	10/Case
Silicone Coated Latex (HCPCS Code A4314)	AS89814	14 Fr.	5 cc	10/Case
	AS89816	16 Fr.	5 cc	10/Case
	AS89818	18 Fr.	5 cc	10/Case
Hydrophilic Coated Latex (HCPCS Code A4314)	AS89814H	14 Fr.	5 cc	10/Case
	AS89816H	16 Fr.	5 cc	10/Case
	AS89818H	18 Fr.	5 cc	10/Case

AMSure® Add-a-Foley Urine Meter Trays

The AMSure Add-a-Foley urine meter kits provide all the supplies necessary for an aseptic insertion. They also offer the caregiver the flexibility to choose the best Foley catheter for the patient and the situation, along with providing a urine meter to obtain the most accurate urine output measurements. **(HCPCS Code A4354)**

Trays include: 400 ml Urine Meter Preconnected with 5 cc Foley Catheter, waterproof and fenestrated drapes, (2) Nitrile gloves, 10% Povidone Iodine solution, (5) Rayon balls, forceps, lubricating jelly, syringe prefilled with 10 cc sterile water, 2 oz. specimen container with a patient label, central supply room wrap, Tyvek header bag, Sterile, not manufactured with rubber latex.



AS89800

ITEM CODE	DESCRIPTION	PRE-FILLED SYRINGE	QTY/UNIT
AS89800	Add-a-Foley Tray 400 mL Urine Meter with 2600 mL Urine Drainage Bag	10 cc	10/Case
AS89830	Add-a-Foley Tray 400 mL Urine Meter with 2600 mL Urine Drainage Bag	30 cc	10/Case



AMSure® Foley Insertion Trays

AMSure Foley Insertion Trays improve simplicity and efficiency of care by providing the basic supplies required to safely insert a Foley catheter. **(HCPCS code A4310)**

- Deep tray basin designed to prevent contamination
- Gloves and povidone iodine swab sticks help maintain sterility and minimized risk of infection
- Sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex

Trays include: 1000 mL collection basin, (3) Povidone-iodine swab sticks or benzalkonium chloride prep pad, vinyl powder-free gloves, waterproof underpad 13 in. x 18 in. (33.02 cm x 45.72 cm), lubrication jelly, fenestrated drape.

SWAB STICK	ITEM CODE	PRE-FILLED SYRINGE	QTY/UNIT
PVP Swab Sticks	AS880	10 cc	10/Case
	AS890	30 cc	10/Case
BZK Swab Sticks	AS880K	10 cc	10/Case
	AS890K	30 cc	10/Case



Prevention of Catheter Associated Urinary Tract Infections


Prevention of Catheter Associated Urinary Tract Infections (CAUTIs) is one of the highest concerns whenever a patient is catheterized for long periods of time. In fact, the longer the period of catheterization, the higher the incidence rate of UTIs. According to the US Centers for Disease Control (CDC), UTIs account for over 30% of infections reported by acute care hospitals.¹ Virtually all of these are caused by introduction of catheters or other instruments into the urinary tract.

The use of a closed system with a sterile fluid path, such as Amsino's AMSure® Pre-Connected Foley Trays, aids in reducing CAUTIs by following the CDC Guidelines² in helping to eliminate breaks in aseptic technique and disconnections of the catheter from the drainage tube.


The AMSure closed Foley catheter trays are not only preconnected at the catheter-tubing junctions, but also include tamper-evident tape at the junction to aid in revealing if the system has been opened.

Furthermore, evidence points to the use of hydrophilic coated catheters as beneficial to patient in comfort, and helping to reduce trauma and infections.^{3,4} Amsino's AMSure Hydrophilic-Coated Foley catheters utilize a proprietary coating method that offer superior lubricity.

See pages 1 and 3 for closed system AMSure Pre-Connected Foley Trays.



Before coating



After coating

Water drop on the surface

AMSure® Urinary Drainage Bags

AMSure urine drainage bags are available in a variety of different configurations to meet the clinicians' needs to best provide patient safety and comfort while helping to minimize infections. All drainage bags include a sheet clamp for the positioning and securing of the drainage tube. The drainage tube has an ID of 11/32" (0.87 cm) for unrestricted flow and are 4 feet (1.2 M) in length, terminating with a standard terraced catheter fitting for a secure and leak-proof connection. Each bag includes a storage pocket for the emptying tube.

(HCPCS code A4358)


- 2,000 cc capacity graduated in 200 cc increments
- Anti-reflux chambers help improve patient safety by restricting urine backflow
- Leak-proof seams ensure fluid containment
- Needle-free sampling port reduces the risk of needle stick injuries
- Air vents help to maintain consistent urine flow, delivering safer patient outcomes
- Sterile fluid pathway, BPA-free, not manufactured with natural rubber latex




Anti-Reflux Drip Device



Anti-Reflux Large Chamber



Anti-Reflux Small Chamber



AMSafe® Needleless Sampling Port, Luer Lock or Luer Slip



Pre-Pierced Blunt Cannula Port



T-Tap Drain Port

ITEM CODE	BAG SHAPE	ANTI-REFLUX TYPE	NEEDLE-FREE SAMPLING PORT	DRAIN TYPE	HANGER TYPE	QTY/UNIT
AS312	Diamond	Drip Device	Pre-Pierced Blunt Cannula	T-Tap	Double-Hook & Rope	20/Case
AS322	Teardrop	Drip Device	Pre-Pierced Blunt Cannula	T-Tap	Double-Hook & Rope	20/Case
AS332	Teardrop	Small Chamber	Pre-Pierced Blunt Cannula	T-Tap	Double-Hook & Rope	20/Case
AS302	Pear	Flutter Valve	Pre-Pierced Blunt Cannula	T-Tap	Single-Hook	20/Case
AS32200	Low-Profile	Large Chamber	Pre-Pierced Blunt Cannula	T-Tap	Swivel Single-Hook & Rope	20/Case
AS32600*	Low-Profile	Large Chamber	Male Luer-Lock or Luer-Slip	T-Tap	Swivel Single-Hook & Rope	20/Case
AS326*	Oval	Drip Device	Pre-Pierced Blunt Cannula	Metal Clamp with Doghouse Storage	Double-Hook & Rope	20/Case

* 100% sterile packed in peel pouch.



AS312 Diamond



AS302 Pear



AS322 Teardrop



AS32200 Low-Profile




AS32600 Low-Profile



AS332 Teardrop



AS326 Oval



AMSure® Urine Meter

AMSure Urine Meter Drain Bags are designed to provide all caregivers the ability to easily and accurately measure urine output. The meter device is pre-attached to a urine drainage bag, for convenient draining of the meter while still obtaining incremental urine measurements. **(HCPCS code A4358)**

- Anti-Reflux Chamber helps improve patient safety by reducing urine backflow
- Twist valve drains at the bottom of the urine meter for easy drainage
- Leak-proof seams ensure fluid containment
- Needle-free sampling port reduces the risk of needle stick injuries
- Air vents help to maintain consistent urine flow delivering safer patient outcomes
- Sterile, not manufactured with natural rubber latex


ITEM CODE	DESCRIPTION	QTY/UNIT
AS328	400 mL Urine Meter with 2600 mL Urine Drainage Bag	10/Case

ACCU-MEASURE® Pediatric Urine Meter

2000 mL Bag with 150 mL Chamber

ACCU-MEASURE Pediatric Urine Meter is specifically designed for pediatric patients. With 1 mL increments, the 150 mL meter chamber's white background provides accurate and precise output measurements. The needle-free sample port reduces the risk of needle stick injuries, helping keep nurses safe while making their job easier. (HCPCS code A4358)

- Double outlet valves – Twist valve for taking sample, T-valve for emptying the bag
- Needle-free sampling port with luer lock access
- Leak-proof seams ensure fluid containment
- Meter bottle capacity 150 mL
- Adjustable bed sheet clamp for securing drainage bag tubing
- Sterile, BPA-free, not manufactured with natural rubber latex



ITEM CODE	DESCRIPTION	QTY/UNIT
AS327100	ACCU-MEASURE Pediatric Urine Meter 2000 mL Bag with 150 mL Chamber	30/Case

Foley Catheters

AMSure® Hydrophilic Latex Foley Catheters

(HCPCS code A4338)

- Hydrophilic coating provides for lubricious surface when contacting moisture to help reduce irritation of the urethra
- Smooth eyes for improved patient comfort
- Symmetrical balloon provides a tight seal at the neck of the urethra
- Sterile, DEHP-free, BPA-free



BALLOON SIZE	ITEM CODE	FOLEY SIZE	QTY/UNIT
5 cc	AS41014H	14 Fr. × 16 in. (40.6 cm)	10/Box
	AS41016H	16 Fr. × 16 in. (40.6 cm)	10/Box
	AS41018H	18 Fr. × 16 in. (40.6 cm)	10/Box
30 cc	AS42014H	14 Fr. × 16 in. (40.6 cm)	10/Box
	AS42016H	16 Fr. × 16 in. (40.6 cm)	10/Box
	AS42018H	18 Fr. × 16 in. (40.6 cm)	10/Box

AMSure® Silicone-Coated Latex Foley Catheters

(HCPCS code A4338)

- Silicone coating helps reduce urethral irritation
- Smooth beveled symmetrical eyes minimize tissue damage
- Large oval eyes aid in maximum drainage
- Firm catheter tip designed to improve ease of insertion
- Symmetrical balloon provides a tight seal at the neck of the urethra
- Individually packaged, sterile, DEHP-free, BPA-free



BALLOON SIZE	ITEM CODE	FOLEY SIZE	QTY/UNIT
3 cc	AS40008	8 Fr. × 10.6 in. (27cm)	10/Box
	AS40010	10 Fr. × 10.6 in. (27cm)	10/Box
5 cc	AS41012	12 Fr. × 16 in. (40.6cm)	10/Box
	AS41014	14 Fr. × 16 in. (40.6cm)	10/Box
	AS41016	16 Fr. × 16 in. (40.6cm)	10/Box
	AS41018	18 Fr. × 16 in. (40.6cm)	10/Box
	AS41020	20 Fr. × 16 in. (40.6cm)	10/Box
	AS41022	22 Fr. × 16 in. (40.6cm)	10/Box
	AS41024	24 Fr. × 16 in. (40.6cm)	10/Box
	AS41026	26 Fr. × 16 in. (40.6cm)	10/Box
	AS41028	28 Fr. × 16 in. (40.6cm)	10/Box
	AS41030	30 Fr. × 16 in. (40.6cm)	10/Box
30 cc	AS42012	12 Fr. × 16 in. (40.6cm)	10/Box
	AS42014	14 Fr. × 16 in. (40.6cm)	10/Box
	AS42016	16 Fr. × 16 in. (40.6cm)	10/Box
	AS42018	18 Fr. × 16 in. (40.6cm)	10/Box
	AS42020	20 Fr. × 16 in. (40.6cm)	10/Box
	AS42022	22 Fr. × 16 in. (40.6cm)	10/Box
	AS42024	24 Fr. × 16 in. (40.6cm)	10/Box
	AS42026	26 Fr. × 16 in. (40.6cm)	10/Box
	AS42028	28 Fr. × 16 in. (40.6cm)	10/Box
	AS42030	30 Fr. × 16 in. (40.6cm)	10/Box

AMSure® 100% Silicone Foley Catheters

(HCPCS code A4344)

- Eliminates latex allergy concerns
- Smooth eyes improve patient comfort
- Reinforced tip provides a stylet for optimized insertion efficacy
- 100% silicone material has non-adherent characteristics and reduces encrustation and calcifications
- Symmetrical balloon provides a tight seal at the neck of the urethra
- Individually packaged, sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex



BALLOON SIZE	ITEM CODE	FOLEY SIZE	QTY/UNIT
5 cc	AS41012S	12 Fr. × 16 in. (40.6cm)	10/Box
	AS41014S	14 Fr. × 16 in. (40.6cm)	10/Box
	AS41016S	16 Fr. × 16 in. (40.6cm)	10/Box
	AS41018S	18 Fr. × 16 in. (40.6cm)	10/Box
	AS41020S	20 Fr. × 16 in. (40.6cm)	10/Box
	AS41022S	22 Fr. × 16 in. (40.6cm)	10/Box
	AS41024S	24 Fr. × 16 in. (40.6cm)	10/Box
	AS41026S	26 Fr. × 16 in. (40.6cm)	10/Box
	AS41028S	28 Fr. × 16 in. (40.6cm)	10/Box
30 cc	AS42012S	12 Fr. × 16 in. (40.6cm)	10/Box
	AS42014S	14 Fr. × 16 in. (40.6cm)	10/Box
	AS42016S	16 Fr. × 16 in. (40.6cm)	10/Box
	AS42018S	18 Fr. × 16 in. (40.6cm)	10/Box
	AS42020S	20 Fr. × 16 in. (40.6cm)	10/Box

Intermittent Urinary Catheters



AS870

AMSure® Urethral Insertion Trays and Kits

AMSure Urethral Trays improve simplicity and efficiency of care by including all components needed to perform short-term aseptic catheterization. (HCPCS code A4353)

- Gloves and povidone iodine swab sticks help maintain aseptic technique and reduce the risk of infection

PVC Catheter Tray

- Staggered eyes help maintain tip strength during insertion and promote efficient drainage
- Short-term catheter use helps decrease infections
- Sterile, BPA-free, not manufactured with natural rubber latex

Red Rubber Catheter Trays

- Smooth tip and quality rubber latex promote patient comfort during insertion
- Clearly marked French size on funnel promotes ease of selection
- Sterile, DEHP-free, BPA-free

Trays include: 1000 mL outer basin tray, 16 in. urethral catheter, waterproof and fenestrated drapes, vinyl powder-free gloves, (3) 10% povidone iodine swab sticks, lubricating jelly, 4 oz. specimen container

ITEM CODE	URETHRAL CATHETER TYPE	QTY/UNIT
AS870	14 Fr PVC*	20/Case
AS87114	14 Fr Red Rubber Latex	20/Case
AS87118	18 Fr Red Rubber Latex	20/Case

*Not manufactured with natural rubber latex.

AMSure® Closed System Urethral Drainage Kits

AMSure Pre-connected Drainage Bag Urethral Catheter Kits improve simplicity and efficiency of care. The catheter is pre-attached to the drainage bag which helps to minimize contamination.

(HCPCS code A4353)

- Gloves and povidone iodine swab sticks help maintain sterility and reduce the risk of infection
- Sterile

Red Rubber Catheter Kits

- Smooth tip and quality rubber latex promote patient comfort during insertion
- Clearly marked French size on funnel promotes ease of selection
- Sterile, DEHP-free, BPA-free

PVC Catheter Kits

- Staggered eyes help maintain tip strength during insertion and promote efficient drainage
- Short-term catheter use helps decrease infections
- Sterile, BPA-free, not manufactured with natural rubber latex



AS89214R

Kit includes: 1200 mL urinary drainage bag, 16 in. (40.64 cm) urethral catheter, CSR wrap, waterproof and fenestrated drapes, (2) Nitrile powder-free gloves, (3) Povidone-iodine swab sticks, lubricating jelly, 2 oz. specimen container with ID label

SWAB STICK	ITEM CODE	FOLEY SIZE	FOLEY MATERIAL	QTY/UNIT
PVP Swab Sticks	AS89214R	14 Fr. × 16 in. (40.6cm)	Red Rubber (Contains Latex)	20/Case
	AS89214	14 Fr. × 16 in. (40.6cm)	PVC Latex- Free	20/Case
BZK Swab Sticks	AS89214RK	14 Fr. × 16 in. (40.6cm)	Red Rubber (Contains Latex)	20/Case
	AS89216RK	16 Fr. × 16 in. (40.6cm)	Red Rubber (Contains Latex)	20/Case

AMSure® Self-Catheterization Kits with R))POLISHED™ Eye Catheters

AMSure Self-Catheterization kits with R-Polish eyelets technology allow the patient, even with limited dexterity, to perform incontinence care with confidence. R-Polish technology uses radio frequency to create ultra-smooth eyelets for optimal patient comfort.

The pre-lubricated, closed system with introducer tip helps to reduce potential for infections. The kits contain all the necessary components needed to perform aseptic self-catheterization. **(HCPCS code A4353)**

- Advanced R-Polished eyelets enhance patient comfort during catheter insertion and removal
- Introducer tip reduces the potential for urinary tract infection (UTI)
- Collection bag with built-in thumb holes allows for ease of use
- Sterile, BPA-free, not manufactured with natural rubber latex

Kit includes: Pre-lubricated 16 in. (40.64 cm) long PVC urethral catheter, 1100 mL collection bag with insertion tip and protector cap, (3) Povidone-iodine swab sticks and benzalkonium chloride prep pad, vinyl powder-free gloves, waterproof underpad 13 in. x 18 in. (33.02 cm x 45.72 cm), waste disposal bag, gauze pad 4 in. x 4 in. (10.16 cm x 10.16 cm)



SWAB STICK	ITEM CODE	CATHETER SIZE	QTY/UNIT
PVP Swab Sticks	AS85008	8 Fr. x 16 in. (40.6cm)	50/Case
	AS85010	10 Fr. x 16 in. (40.6cm)	50/Case
	AS85012	12 Fr. x 16 in. (40.6cm)	50/Case
	AS85014	14 Fr. x 16 in. (40.6cm)	50/Case
	AS85016	16 Fr. x 16 in. (40.6cm)	50/Case
BZK Swab Sticks	AS85008K	8 Fr. x 16 in. (40.6cm)	50/Case
	AS85010K	10 Fr. x 16 in. (40.6cm)	50/Case
	AS85012K	12 Fr. x 16 in. (40.6cm)	50/Case
	AS85014K	14 Fr. x 16 in. (40.6cm)	50/Case
	AS85016K	16 Fr. x 16 in. (40.6cm)	50/Case

AMSure® PVC Intermittent Urethral Catheter with R))POLISHED™ Eyes

AMSure PVC Intermittent Catheters with R-Polish technology are designed for short-term use. R-Polish technology uses radio frequency to create ultra-smooth eyelets for optimal patient comfort. The intermittent catheters are available in various French sizes as well as female, pediatric and male lengths. (HCPCS code A4351)

- Advanced R-Polished eyelets enhance patient comfort during catheter insertion and removal
- Staggered eyes help maintain tip strength during insertion and promote efficient drainage
- Designed to minimize trauma and reduce risk of infection
- Sterile, BPA-free, not manufactured with natural rubber latex



FEATURES	ITEM CODE	DESCRIPTION	QTY/UNIT
Female	AS960610	10 Fr, 6 in. (15.24 cm)	50/Case
	AS960612	12 Fr, 6 in. (15.24 cm)	50/Case
	AS960614	14 Fr, 6 in. (15.24 cm)	50/Case
	AS960616	16 Fr, 6 in. (15.24 cm)	50/Case
Pediatric	AS961010	10 Fr, 10 in. (25.40 cm)	50/Case
Male	AS961610	10 Fr, 16 in. (40.64 cm)	50/Case
	AS961612	12 Fr, 16 in. (40.64 cm)	50/Case
	AS961614	14 Fr, 16 in. (40.64 cm)	50/Case
	AS961616	16 Fr, 16 in. (40.64 cm)	50/Case
Coudé Tip	AS861610C	10 Fr, 16" (40.65 cm)	30/Case
	AS861612C	12 Fr, 16" (40.65 cm)	30/Case
	AS861614C	14 Fr, 16" (40.65 cm)	30/Case
	AS861616C	16 Fr, 16" (40.65 cm)	30/Case
	AS861618C	18 Fr, 16" (40.65 cm)	30/Case

AMSure® PVC Intermittent Urethral Catheter

AMSure PVC Intermittent Urethral Catheters provide a smooth surface firm enough for insertion yet soft enough to be slightly pliable. The urethral catheters are designed for intermittent use and available in various French sizes as well as female, pediatric and male lengths. **(HCPCS code A4351)**

- Staggered eyes help maintain tip strength during insertion and promote efficient drainage
- Short-term catheter use helps decrease infections
- Clear shaft provides easy viewing as you measure
- PVC material eliminates latex allergy concerns
- Individually packaged, sterile, BPA-free, not manufactured with natural rubber latex

FEATURES	ITEM CODE	DESCRIPTION	QTY/UNIT
Female	AS860608	8 Fr, 6 in. (15.24 cm)	50/Case
	AS860610	10 Fr, 6 in. (15.24 cm)	50/Case
	AS860612	12 Fr, 6 in. (15.24 cm)	50/Case
	AS860614	14 Fr, 6 in. (15.24 cm)	50/Case
	AS860616	16 Fr, 6 in. (15.24 cm)	50/Case
Pediatric	AS861006	6 Fr, 10 in. (25.40 cm)	50/Case
	AS861008	8 Fr, 10 in. (25.40 cm)	50/Case
	AS861010	10 Fr, 10 in. (25.40 cm)	50/Case
Male	AS861608	8 Fr, 16 in. (40.64 cm)	50/Case
	AS861610	10 Fr, 16 in. (40.64 cm)	50/Case
	AS861612	12 Fr, 16 in. (40.64 cm)	50/Case
	AS861614	14 Fr, 16 in. (40.64 cm)	50/Case
	AS861616	16 Fr, 16 in. (40.64 cm)	50/Case
	AS861618	18 Fr, 16 in. (40.64 cm)	50/Case





AMSure® Red Rubber Urethral Catheter

AMSure Red Rubber Urethral Catheters are pliable and designed for intermittent use. The urethral catheters are available in various French sizes as well as female and male lengths.

(HCPCS code A4351)

- Smooth tip and quality rubber latex promote patient comfort during insertion
- Clearly marked French size on funnel promotes ease of selection
- Individually packaged, sterile, DEHP-free


ITEM CODE	DESCRIPTION	QTY/UNIT
AS44008	8 Fr, 11 in. (27.94 cm)	100/Case
AS44010	10 Fr, 11 in. (27.94 cm)	100/Case
AS44012	12 Fr, 16 in. (40.64 cm)	100/Case
AS44014	14 Fr, 16 in. (40.64 cm)	100/Case
AS44016	16 Fr, 16 in. (40.64 cm)	100/Case
AS44018	18 Fr, 16 in. (40.64 cm)	100/Case
AS44020	20 Fr, 16 in. (40.64 cm)	100/Case
AS44022	22 Fr, 16 in. (40.64 cm)	100/Case
AS44024	24 Fr, 16 in. (40.64 cm)	100/Case
AS44028	28 Fr, 16 in. (40.64 cm)	100/Case

Irrigation Sets

AMSure® Cystoscopy & TUR Irrigation Sets

The AMSure Cystoscopy and TUR Irrigation sets are used for continuous bladder irrigation after prostate or bladder surgeries. These are available with one, two, or four leads to meet a variety of needs. (HCPCS code A4351)

- Robert Clamp(s)
- Drip Chamber(s)
- Sterile, not manufactured with natural rubber latex



ITEM CODE	TYPE	QTY/UNIT
ABL400	Single Lead	20/Case



ITEM CODE	TYPE	QTY/UNIT
ABL402	Two Leads	20/Case



ITEM CODE	TYPE	QTY/UNIT
ABL404	Four Leads	20/Case

Ancillary & Patient Care Products

Urinary Leg Bags

(HCPCS code A4358)

- Adjustable, soft, breathable pre-attached elastic leg straps for patient comfort
- Anti-reflux flutter valve improves patient safety by reducing urine backflow
- Leak-proof seams ensure fluid containment
- Sterile fluid pathway
- BPA-free, not manufactured with natural rubber latex



AS307N

SIZE	ITEM CODE	DRAIN PORT STYLE	18 in. (45.72 cm) EXTENSION TUBE WITH ADAPTER & PROTECTIVE CAP	QTY/UNIT
600 mL	AS306N	Push-Pull		48/Case
	AS308N	Twist-Turn		48/Case
	LBE306N	Push-Pull	•	48/Case
	LBE308N	Twist-Turn	•	48/Case
900 mL	AS307N	Push-Pull		48/Case
	AS309N	Twist-Turn		48/Case
	LBE307N	Push-Pull	•	48/Case
	LBE309N	Twist-Turn	•	48/Case

Extension Tube for Leg Bags

(HCPCS code A4331)

- Kink-resistant tubing safeguards against fluid buildup
- Flexible tubing permits leg bag to be positioned on the thigh or calf
- Sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex



ITEM CODE	DESCRIPTION	QTY/UNIT
AS311	18 in. (45.72 cm) Extension Tube	50/Case



Pediatric Urine Collector

(HCPCS code A4328)

- Adhesive surface folds to form a spill-proof specimen transport container
- Hypoallergenic feature makes it safe to use with most pediatric patients
- Visible graduation marks designed for ease of use
- Single-use disposable
- Sterile, DEHP-free, BPA-free, not manufactured with natural rubber latex

ITEM CODE	DESCRIPTION	QTY/UNIT
AS409	200 mL Pediatric Urine Bag	50/Box



AMSure® Irrigation Trays and Kits

AMSure Irrigation Kits and Trays are designed for simplicity and efficiency of care. The trays are available sterile or non-sterile to suit your clinical need.

(HCPCS code A4320)

Bulb Syringe

- Flexible soft bulb promotes one-handed administration
- Designed to stand on the bulb end to help minimize contamination
- Patient label provides convenient patient documentation
- DEHP-free, BPA-free, not manufactured with natural rubber latex
- Includes syringe tip protector

AS121



Thumb Control Ring Syringe

- Thumb control ring top allows for maximum control and one-handed administration
- Precision crafted O-ring plunger ensures a positive seal
- Re-sealable pole bag designed for safe and convenient storage
- Patient label provides convenient patient documentation
- DEHP-free, BPA-free, not manufactured with natural rubber latex
- Includes syringe tip adapter and protector cap

AS126



Non-sterile irrigation kits include: 500 mL graduated container, patient ID label and tip protector.

Sterile irrigation trays include: 1000 mL outer tray, 500 mL graduated container, waterproof underpad, alcohol prep pad, patient ID label and tip protector.

Non-sterile Irrigation Trays

ITEM CODE	THUMB CONTROL RING SYRINGE	BULB SYRINGE	RESEALABLE POLE BAG	QTY/UNIT
AS121		•		30/Case
AS126	•			30/Case
AS127	•		•	30/Case

Sterile Irrigation Trays

ITEM CODE	THUMB CONTROL RING SYRINGE	BULB SYRINGE	QTY/UNIT
AS130		•	20/Case
AS136	•		20/Case

AS136



Syringes

(HCPCS code A4322)

Bulb Syringe

- Flexible soft bulb promotes one-handed administration
- Designed to stand on the bulb end to help minimize contamination
- DEHP-free, BPA-free, not manufactured with natural rubber latex
- Includes syringe tip protector

Piston Thumb Control Ring Syringe

- Piston thumb control ring top allows for maximum control and one-handed administration
- Precision crafted O-ring plunger ensures a positive seal
- Re-sealable pole bag designed for safe and convenient storage
- DEHP-free, BPA-free, not manufactured with natural rubber latex
- Includes syringe luer tip adapter and protector cap

Flat Top Syringe

- Visible dosing marks designed for ease of use
- Precision crafted O-ring plunger ensures a positive seal
- Visible dosing marks designed for ease of use
- DEHP-free, BPA-free, not manufactured with natural rubber latex
- Available in sterile or non-sterile packaging
- Includes syringe luer tip adapter and protector cap



TYPE	ITEM CODE	PACKAGING	STERILITY	TIP PROTECTOR CAP	SMALL TUBE ADAPTER	QTY/UNIT
Bulb	AS011P	Form Fill & Seal	Sterile (EO)	•	•	50/Case
	AS011	Poly Pouch	Sterile (EO)	•		50/Case
	AS010	Bulk	Non-Sterile			50/Case
Piston Thumb Control Ring	AS015	Poly Pouch	Sterile (EO)	•	•	50/Case
	AS016	Resealable Hanger Bag	Non-Sterile	•	•	30/Case
Piston Flat Top	AS115	Poly Pouch	Sterile (EO)	•	•	50/Case
	AS116	Resealable Pole Bag	Non-Sterile	•	•	30/Case

AMSure® Specimen Containers

AMSure 4 oz., leak-resistant specimen containers are designed for the collection and safe transport of patient samples. (HCPCS code P9615)

- DEHP-free, BPA-free, not manufactured with natural rubber latex



AS340

ITEM CODE	PACKAGING	STERILITY	LABELS	QTY/UNIT
AS340	Poly Pouch	Interior	Patient Information	100/Case
AS341	Bulk	Non-Sterile	Tamper Evident	100/Case
AS343	Bulk	Non-Sterile	None	500/Case
AS345	Form-Fill & Seal	Fully Sterile	None	100/Case

Index

ITEM CODE	HCPCS CODE	PAGE NUMBER
ABL400	A4351	20
ABL402	A4351	20
ABL404	A4351	20
AS010	AS4322	24
AS011	AS4322	24
AS011P	AS4322	24
AS015	AS4322	24
AS016	AS4322	24
AS115	AS4322	24
AS116	AS4322	24
AS121	AS4320	23
AS126	AS4320	23
AS127	AS4320	23
AS130	AS4320	23
AS136	AS4320	23
AS302	A4358	7
AS306N	A4358	21
AS307N	A4358	21
AS308N	A4358	21
AS309N	A4358	21
AS311	AS4331	21
AS312	A4358	7
AS322	A4358	7
AS32200	A4358	7
AS326	A4358	7
AS32600	A4358	7
AS327100	A4358	10
AS328	A4358	8
AS332	A4358	7
AS340	P9615	25
AS341	P9615	25
AS343	P9615	25
AS345	P9615	25
AS40008	A4338	12
AS40010	A4338	12
AS409	AS4328	22
AS41012S	A4344	13
AS41012	A4338	12
AS41014	A4338	12
AS41014S	A4344	13
AS41014H	A4338	11

ITEM CODE	HCPCS CODE	PAGE NUMBER
AS41016	A4338	12
AS41016S	A4344	13
AS41016H	A4338	11
AS41018	A4338	12
AS41018S	A4344	13
AS41018H	A4338	11
AS41020	A4338	12
AS41020S	A4344	13
AS41022	A4338	12
AS41022S	A4344	13
AS41024	A4338	12
AS41024S	A4344	13
AS41026	A4338	12
AS41026S	A4344	13
AS41028	A4338	12
AS41028S	A4344	13
AS41030	A4338	12
AS42012	A4338	12
AS42012S	A4344	13
AS42014	A4338	12
AS42014S	A4344	13
AS42014H	A4338	11
AS42016	A4338	12
AS42016S	A4344	13
AS42016H	A4338	11
AS42018	A4338	12
AS42018S	A4344	13
AS42018H	A4338	11
AS42020	A4338	12
AS42020S	A4344	13
AS42022	A4338	12
AS42024	A4338	12
AS42026	A4338	12
AS42028	A4338	12
AS42030	A4338	12
AS44008	A4351	19
AS44010	A4351	19
AS44012	A4351	19
AS44014	A4351	19
AS44016	A4351	19
AS44018	A4351	19

Index

ITEM CODE	HCPCS CODE	PAGE NUMBER
AS44020	A4351	19
AS44022	A4351	19
AS44024	A4351	19
AS44028	A4351	19
AS85008	A4353	15
AS85008K	A4353	16
AS85010	A4353	16
AS85010K	A4353	16
AS85012	A4353	16
AS85012K	A4353	16
AS85014	A4353	16
AS85014K	A4353	16
AS85016	A4353	16
AS85016K	A4353	16
AS860608	A4351	18
AS860610	A4351	18
AS860612	A4351	18
AS860614	A4351	18
AS860616	A4351	18
AS861006	A4351	18
AS861008	A4351	18
AS861010	A4351	18
AS861608	A4351	18
AS861610	A4351	18
AS861610C	A4351	17
AS861612	A4351	18
AS861612C	A4351	17
AS861614	A4351	18
AS861614C	A4351	17
AS861616	A4351	18
AS861616C	A4351	17
AS861618	A4351	18
AS861618C	A4351	17
AS870	A4353	14
AS87114	A4353	14
AS87118	A4353	14
AS880	A4310	5
AS880K	A4310	5
AS890	A4310	5
AS890K	A4310	5
AS89214	A4353	15

ITEM CODE	HCPCS CODE	PAGE NUMBER
AS89214R	A4353	15
AS89214RK	A4353	15
AS89216RK	A4353	15
AS89300	A4354	2
AS89314	A4314	1
AS89314H	A4314	1
AS89314S	A4315	1
AS89316	A4314	1
AS89316H	A4314	1
AS89316S	A4315	1
AS89318	A4314	1
AS89318H	A4314	1
AS89318S	A4315	1
AS89330	A4354	2
AS89800	A4354	4
AS89814	A4314	3
AS89814H	A4314	3
AS89814S	A4315	3
AS89816	A4314	3
AS89816H	A4314	3
AS89816S	A4315	3
AS89818	A4314	3
AS89818H	A4314	3
AS89818S	A4315	3
AS89830	A4354	4
AS960610	A4351	17
AS960612	A4351	17
AS960614	A4351	17
AS960616	A4351	17
AS961010	A4351	17
AS961610	A4351	17
AS961612	A4351	17
AS961614	A4351	17
AS961616	A4351	17
LBE306N	A4358	21
LBE307N	A4358	21
LBE308N	A4358	21
LBE309N	A4358	21

Sources

1. Centers for Disease Control and Prevention, <https://www.cdc.gov/infectioncontrol/guidelines/cauti/background.html>. Published November 5, 2015. Accessed September 30, 2022.
2. Gould, MD, MSCR, C., Umscheid, MD, MSCE, C., Agarwal, MD, MPH, R., Kuntz, MSW, MSLIS, G., Pegues, MD, D., & Healthcare Infection Control Practices Advisory Committee (HICPAC). (2019, June 6). GUIDELINE FOR PREVENTION OF CATHETER ASSOCIATED URINARY TRACT INFECTIONS 2009. [Www.Cdc.Gov/Infectioncontrol/Pdf/Guidelines/Cauti-Guidelines-H.Pdf](https://www.cdc.gov/infectioncontrol/Pdf/Guidelines/Cauti-Guidelines-H.Pdf). Retrieved September 30, 2022, from <https://www.cdc.gov/infectioncontrol/guidelines/cauti/>
3. De Ridder, D. J., Everaert K, Fernandez LG, et al. Intermittent catheterisation with hydrophilic-coated catheters (SpeediCath) reduces the risk of clinical urinary tract infection in spinal cord injured patients: A prospective randomised parallel comparative trial. *EurUrol.* 2005;48(6):991-995.
4. Stensballe J, Looms D, Nielsen PN, Tvede M. Hydrophilic-coated catheters for intermittent catheterisation reduce urethral micro trauma: A prospective, randomised, participant-blinded, crossover study of three different types of catheters. *EurUrol.* 2005;48(6):978-983.



Advancing Healthcare Worldwide[®]

Amsino Medical Group Co. Ltd.

708 Corporate Center Drive
Pomona, CA 91768
USA

P: +1 800.632.6746

F: +1 909.626.3888

info@amsino.com

www.amsino.com