




**SINGLE USE  
FLEXIBLE  
CYSTOSCOPE**



Shanghai AnQing Medical Instrument Co., Ltd.  
3&4 Floor, No.2 Building, 366 Huiqing Rd,  
East Zhangjiang High-Tech Park,  
201201 Shanghai, China


Tel: +86 21 61117375  
Fax: +86 21 61117374  
Email: [global@innovexmed.com.cn](mailto:global@innovexmed.com.cn)

**SINGLE USE STILL  
NEEDS PERFORMANCE**

# BRIGHT & CLEAR

## HIGH IMAGE QUALITY


INNOVEX Single-Use Flexible Cystoscope provides every patient a clean scope with high quality view - every time. Eliminate the cost and risk associated with reprocessing and cross-contamination.



Metal snake bone structure  
Excellent maneuverability  
Working length 380mm

Dual deflection  
210°up / 180°down

Latest CMOS technology  
Superior image  
Field of view 110°±10%  
Depth of field: 5mm-100mm



## CLEAR & SHARP IMAGE

DMR-UB1-013-01 A/0