

SINGLE USE FLEXIBLE URETERORENOSCOPE



SHANGHAI ANQING MEDICAL INSTRUMENT CO, LTD.
3&4 Floor, No.2 Building, 366 Huiqing Rd,
East Zhangjiang High-Tech Park,
201201 Shanghai, China


Tel: +86 21 61117375
Fax: +86 21 61117374
Email: global@innovexmed.com.cn

**SINGLE USE STILL
NEED PERFORMANCE**

BRIGHT & CLEAR

HIGH IMAGE QUALITY

INNOVEX Single-Use Flexible Ureterorenoscope provides every patient a clean scope with high quality view – every time. Eliminate the cost and risk associated with reprocessing and cross-contamination.



US 31E-12-EU
US 31E-12-US



CLEAR & SHARP IMAGE

Dual LED illumination
Latest CMOS technology
Superior image
Field of view 110° ± 10%
Depth of field: 3mm-100mm



US 31D-12-EU
US 31D-12-US